



# New Developments Triskem - Part 2

UGM - 21/09/2018

Jesus College - Cambridge (UK)

S. Happel

Radiopharmacy  
and  
Nuclear Medicine

Environment and  
Bioassay

Geochemistry  
and  
Metals Separation

Decommissioning



# New Resins - domains and applications

- Analytical

- Radiochemistry

- Environmental monitoring, bioassay, waste monitoring, decommissioning
- Actinides, fission and activation products, NORM, methods/resins for DTMs, rapid methods,...
  - » **TK100/1, TK200, TBP, CL Resin, TK201/2, TK300, TK-TcScint, Extractive discs,...**




- Mass spectrometry

- Isotope ratio determination (universities, petrol industry,...)
- Sr, Pb, U, actinides, Cu, Sn...
  - » Dating of geological samples : **TK200, TK300, TK400,...**
  - » Food provenancing : **TK100,...**
  - » Nuclear forensics: **TK200, TK400,...**
  - » Biomedical (Cu),... => **CU Resin,...**



- Decommissioning/decontamination
  - Treatment of effluents / liquid wastes / environmental waters
    - Removal of radioactive contaminants & heavy metals (Cs, Sr, Ra, I...)
  - Inorganic compounds embedded into PAN matrix
    - **CS Resins, more under development**
  - Polymer-based => **TK202,...**
- Hydrometallurgy
  - Recovery of critical metals
  - Recycling
    - **Mainly functionalized polymers => under development**



- Radiopharmacy/Nuclear Medicine

Radiopharmacy  
and  
Nuclear Medicine

- Radionuclide production

- Cooperation with cyclotrons & reactors (NL, RN producers,...)

- Separation of radionuclides from irradiated targets

- » Diagnostics: Zr-89, Cu-64, Ga-68, Ge-68, Ti-44/5, Tc-99m...


- **ZR Resin, CU Resin, TK200 Resin, TK400, TK201, TK202,...**

- » Therapy: alpha emitters, Lu-177, Cu-67, Sn-117m,...


- **TK400, TK200, TBP Resin, CU Resin,...**

- Radiopharmacy/nuclear medicine
  - Purification of generator eluates => under development
  - Decontamination of contaminated effluents => **CL Resin**,...
  - Quality control
    - Cartridge based methods
    - **DGA sheets** (functionalized TLC, Ra-223, Ga-68, Pb-212,....)


Radiopharmacy  
and  
Nuclear Medicine





Radiochromatogram measured immediately after separation. Low abundant radiations of  $^{227}\text{Ac}$  were not detected.



Radiochromatogram measured one hour after separation. Decay and ingrowth of  $^{211}\text{Pb}$  is clearly visible.




- Original scope: Hydroxamate based resin
  - Standard for Zr separation from Y targets
  - Ready to use
  - Facile Zr elution (avoid 1M oxalic acid)
- Zr-89 production via (p,n) reaction from natural Y (Y-89) targets
  - High Zr/Y selectivity necessary
- Fast extraction kinetics: > 99% after 1 min
- Also used for other separations possible (Ti/Sc, Ge/Ga, Ga/Zn)
- Applications in analytical applications (Zr-93, Nb,...)



- No selectivity for Y, Sc
- High Ge/Ga selectivity at elevated HCl
- High selectivity for Zr, Ti, Nb, W over wide HCl concentration range

- No selectivity for alkaline and earth alkalines
- Lanthanides are not retained
- Strong Fe retention (dip at 2 – 3M HCl)




- High selectivity for Zr, Ti, Nb, W over wide HNO<sub>3</sub> concentration range
  - Loss of selectivity at 6M HNO<sub>3</sub>  
=> Resin shows colour change
- No selectivity for Y, Sc, lanthanides, earth alkalines, most transition metals,...
- High Ge/Ga selectivity at 3M HNO<sub>3</sub>






# Zr separation incl. Fe removal

## Zr separation on 1 mL ZR Resin column




- Load from 2 – 6M HCl
- Rinsing described by Holland can be used
- No activation with acetonitrile
- Quantitative Zr elution in 1.5 - 2 mL 0.05M oxalic acid
- Clean Fe removal





- Ge separation from Ga, GaNi or GaCo for Ge-68 production
- Loading from HNO<sub>3</sub>, HCl or H<sub>2</sub>SO<sub>4</sub>
  - HNO<sub>3</sub> & H<sub>2</sub>SO<sub>4</sub> preferred
  - => GeCl<sub>4</sub> volatile
- Fe removal with 2M HCl when loading from HNO<sub>3</sub>

Ti/Sc separation, ZR Resin



- In-house tested up to 50 mg Sc (100 mg ZR Resin)
- Ti retained from 10M HCl, Sc not retained
- Ti elution with 0.1M citric, 0.2M oxalic acid, 0.1M H<sub>2</sub>O<sub>2</sub>



- 4g irradiated Sc
- 65.2 MBq Ti-44
- 5 mL Zr Resin
- Ti-44 yield >95%
- $D_f(\text{Sc}): 10^5$

Fig. 3. HCl concentration dependency of  $K_d$  for  $^{44}\text{Ti}/^{46}\text{Sc}$  on ZR hydroxamate resin. Fig. 5.  $^{44}\text{Ti}/^{46}\text{Sc}$  elution profile using ZR hydroxamate resin with a load of 4 g of scandium.

## Use of ZR Resin as support in Ti-44/Sc-44 generators





Fig. 1. Schematic concept of a forward/reverse flow radionuclide generator.

- Direct (1 mL ZR Resin) and reverse elution (2 mL ZR Resin)
- 65 column volumes tested up until publication
- High Sc yields, max. Ti-44 breakthrough:  $4.1 \times 10^{-4}\%$
- Obtained Sc gave labelling yields > 94%
- Generator has just been set-up at BNL


# Ga-67/8 separation from Zn targets

Separation Zr resin - L: 10M HCl - E: 1.5M HCl



EANM '17: Presentation Gagnon et al.

Separation TK200 resin









- Loading from:
  - 0.1M HNO<sub>3</sub> (liquid targets)
  - >6M HCl (solid targets)
- Clean Ga separation on ZR Resin
- Elution with 1.5M HCl
  - Too high for labeling/injection

- Ga conversion step on TK200
- TK200 load from 1,5M HCl, elution in 2 – 3 BV water


# D<sub>w</sub> values TK200 - HCl

=> product sheet on website




# Other examples for separations on TK200


- Pt separation from Ir



- Zn separation from Cu



- Sc from Ti or Ca



More applications under development

- Oxime based resin (very hydrophobic, preferably use as cartridge)
  - Column packing from 20% EtOH
- High selectivity for Cu
- Load from pH >2, elution in high acid
  - Not ideal for solid Ni targets
    - TK201 under evaluation, on its own and combined with CU Resin
  - Very suitable for liquid targets and solid Zn targets (=> Cu-67)
  - High labelling yields
- Typical applications:
  - Cu isotope production from Ni or Zn targets
  - Cu removal (Ni-64 purification)
  - Cu isotope ratio determination by mass spectrometry (=> ESI)
    - Geochemistry, Archaeometry, Biomedical (Wilson's Disease)



## Purification of <sup>67</sup>Cu and Recovery of its Irradiated Zn Target

A.J. DeGraffenreid<sup>a</sup>, R. Nidzyn<sup>a</sup>, B. Jenkins<sup>a</sup>, D.E. Wycoff<sup>b</sup>, T.E. Phelps<sup>b</sup>, A. Goldberg<sup>a</sup>, D.G. Medvedev<sup>a</sup>, S.S. Jurisson<sup>b</sup>, C.S. Cutler<sup>a</sup>

<sup>a</sup>Brookhaven National Laboratory, C-AD/MIRP—Upton, NY (USA)

<sup>b</sup>University of Missouri, Department of Chemistry—Columbia, MO (USA)

Poster  
presented at  
ISRS 2017

### Cu Resin

- Dissolve irradiated (86.4  $\mu$ A for 0.79 d) <sup>nat</sup>Zn target (13.70 g) in trace metal grade conc. HCl
- Evaporated to dryness and reconstitute into milliQ water (91.5 mL) to give a 312 mg ZnCl<sub>2</sub>/mL solution
- Ensure pH  $\geq$  2 with pH paper
- Added ~60.6 mL (18.9 g ZnCl<sub>2</sub>) to Cu Resin column (2.40 g)
- Quantitative transfer with 2 x 20 mL H<sub>2</sub>O  $\rightarrow$  into load bottle
- Rinsed with pH 2 HCl (80 mL)
- Recover radiocopper with 2 x 20 mL 6 M HCl
- Evaporate 6 M HCl to dryness



### Cu Resin

#### Recovery (%)

Nuclide	EOB Activity (mCi $\pm$ 1 $\sigma$ )	Recovery (%)			
		Load w/ Quant. Transfer	pH 2 HCl Rinse	Acid #1	Acid #2
<sup>64</sup> Cu	4700 $\pm$ 200	ND	ND	102	ND
<sup>65</sup> Zn	41.0 $\pm$ 0.8	103	ND	0.04	ND
<sup>58</sup> Co	63 $\pm$ 1	104	0.04	0.1	0.01


- Produced 143 mCi <sup>67</sup>Cu
- Quantitative recovery of radiocopper
- 99.5% radionuclidic purity—single column
- ICP-OES: 132.9  $\mu$ g Cu and 1.3 mg Zn
  - Anion exchange column still needed to remove trace Zn
- Specific activity <sup>67</sup>Cu at EOB: 1.07 mCi/ $\mu$ g

### Cu Resin


Robust separation that could shorten the overall processing time to separate co-produced radionuclides and large quantities of Zn from radiocopper  
Cation and anion exchange columns still needed to suitably purify radiocopper

# Cu separation- conversion on TK201




- Cu conversion from high HCl to more suitable conditions usually done using AIX
- Problem: shrinking/swelling of AIX, insufficient trace Zn removal
- TK201 to be used as replacement?





- Removal of trace Ni during load/rinse
- Zn remains retained during Cu elution
- Solid Ni targets: improvement of Ga/Co removal or combination with CU Resin



- QC of radionuclides and generator eluents (p.ex. Ra-223, Ac-225/Bi-213, Pb-212, Ge-68/Ga-68 ...) with TLC scanner or radiometer/LSC after cutting
  - Therapy: alpha emitters
  - Diagnostics e.g. with generator produced Ga-68
- More types of sheets under development (selectivities, geometry)
- 2D TLC under development (with Subatech) => use in decommissioning
  - Quadratic sheet, two runs (90° turn in-between) with different acids => 2D pattern
  - Measurement e.g. with Ai4r Beaver system (high res  $\alpha/\beta$  discrimination counting)



- Long chained alcohol
- Main application: Pa separation (Pa-230/1 determination by MS)
  - NPL (no selectivity for actinides, Ac, Ra, Pb,...=> Pa-230 purif.)
  - Also retains Mo, Nb, Po, Ga,...
    - Ga-68 from solid Zn targets under testing
    - Nb separation from Zr (Nb-90)
    - Mo separation for Mo isotope ratio determination
  - At separation from Bi targets to replace dry distillation?



## Chromatographic separation of the theranostic radionuclide $^{111}\text{Ag}$ from a proton irradiated thorium matrix

Tara Mastren<sup>a</sup>, Valery Radchenko<sup>a,1</sup>, Jonathan W. Engle<sup>a,2</sup>, John W. Weidner<sup>a</sup>, Allison Owens<sup>b</sup>, Lance E. Wyant<sup>b</sup>, Roy Copping<sup>b</sup>, Mark Brugh<sup>a</sup>, F. Meiring Nortier<sup>a</sup>, Eva R. Birnbaum<sup>a</sup>, Kevin D. John<sup>a</sup>, Michael E. Fassbender<sup>a,\*</sup>


<sup>a</sup> Chemistry Division, Los Alamos National Laboratory, P.O. Box 1663, Los Alamos, NM 87545, USA

<sup>b</sup> Nuclear Security and Isotope Technology Division, Oak Ridge National Laboratory, Oak Ridge, TN 37831, USA

### HIGHLIGHTS

- Chromatographic recovery of medical isotope  $^{111}\text{Ag}$  from proton irradiated thorium targets.
- First-time measured equilibrium distribution coefficients for silver and ruthenium on CL resin.
- $^{232}\text{Th}$  (p, fission) cross-section data for the formation of  $^{111}\text{Ag}$  and  $^{110m}\text{Ag}$ .

### GRAPHICAL ABSTRACT



- On-going tests for Au, Pd, Hg
- Use in DGT



Applied Radiation and Isotopes

Volume 86, April 2014, Pages 90-96



## Accurate determination of $^{129}\text{I}$ concentrations and $^{129}\text{I}/^{137}\text{Cs}$ ratios in spent nuclear resins by Accelerator Mass Spectrometry

Emmanuelle Nottoli<sup>a,\*</sup>, Philippe Bienvenu<sup>b</sup>, Alexandre Labet<sup>a</sup>, Didier Bourlés<sup>c</sup>, Maurice Arnold<sup>c</sup>, Maité Bertaux<sup>a</sup>

Show more

<https://doi.org/10.1016/j.apradiso.2014.01.010>

Get rights and content

Anal. Chem., 2018, <https://pubs.acs.org/doi/10.1021/acs.analchem.8b01380>

## Separation of protactinium employing sulfur-based extraction chromatographic resins

Tara Mastren<sup>†</sup>, Benjamin W. Stein<sup>†</sup>, T. Gannon Parker<sup>†</sup>, Valery Radchenko<sup>†#</sup>, Roy Copping<sup>‡</sup>, Allison Owens<sup>‡</sup>, Lance E. Wyant<sup>‡</sup>, Mark Brugh<sup>†</sup>, Stosh A. Kozimor<sup>†</sup>, F. Meiring Nortier<sup>†</sup>, Eva R. Birnbaum<sup>†</sup>, Kevin D. John<sup>†</sup>, Michael E. Fassbender<sup>†\*</sup>

<sup>†</sup>Chemistry Division, Los Alamos National Laboratory, P.O. Box 1663, Los Alamos, NM 87545, USA

<sup>‡</sup>Nuclear Security and Isotope Technology Division, Oak Ridge National Laboratory, Oak Ridge, TN 37831, USA

<sup>#</sup>Current Address: Life Sciences Division, TRIUMF, 4004 Wesbrook Mall, Vancouver, BC, V6T2A3, Canada

- Decamp et al.: Iodine removal from elevated sample volumes at high flow rates<sup>§</sup>
  - Issues with rad. waste storage
    - Storage as liquid waste challenging
    - Preferably stored as solid waste
  - > 10 L radioactive process effluent (1M HNO<sub>3</sub>),
  - Flow rate up to 180 mL/min
  - Mixed-bed columns
    - 3g CL resin (plus 4g XAD-4 resin), iodine uptake: 85 – 95%

<sup>§</sup> C. Decamp (IRE), S. Happel: Utilization of a mixed-bed column for the removal of iodine from radioactive process waste solutions, Journal of Radioanalytical and Nuclear Chemistry, online April 2013, DOI: 10.1007/s10967-013-2503-1

# Some on-going projects

- C.L.I.P.S. 2020
- Ra Resin(s)
- TK300 Resin for Cs/Rb separation
- Online concentration/separation
- Extractive discs
- Resins based on scintillating beads
  - Uni Barcelona
- Improved Ni Resin
- Li Resin
- Lanthanide separations (e.g. Lu-177 separation from Yb) => new resins
- SE Resin: e.g. Se-72/As-72 generator (piaszelenol chemistry)
- Tc-99m separation from Mo (TK201/2)
- DGT (radioelements and heavy metals)
  - Passive sampling/'bio-availability'
- Ga purification (generator) and QC
- At-211 separation (e.g. TK400)
- Ac-225 separation e.g. from Ra-226
- Removal of radionuclides from effluents
- Functionalised polymers/silicates,...
- ...



Thank you for your attention!

We hope you have enjoyed the UGM  
Thank you for your participation!



Radiopharmacy  
and  
Nuclear Medicine

Environment and  
Bioassay

Geochemistry  
and  
Metals Separation

Decommissioning

