



TrisKem International

An overview over some new extraction chromatographic resins and their application in radiopharmacy

Steffen Happel

04/06/19



TrisKem International



- Based in Rennes (France)
- Independent company since 02/07
 - Formerly part of Eichrom Europe
 - ISO 9001 since 2007
- Staff : 23
- R&D and TechSupport group:
 - 3 RadChem PhD, 2 OrgChem PhD, 1 Engineer and 3 Technicians
- R&D: Development of new resins, techniques and applications
- Several domains:



Radiopharmacy
and
Nuclear Medicine

Environment and
Bioassay

Geochemistry
and
Metals Separation

Decommissioning



- Two R&D labs:

- Synthesis Lab (new resins and extractants)

- Incl. grafted resins (silica or polymers), macrocycles,...

- Application Lab

- Preparation of extraction chromatographic resins
- Resin characterisation and method development



- Equipment:

- ICP-MS, IC, TOC, TGA, IR, automatic desiccator, benchtop NMR (43 MHz), surface area and pore size volume analyser, size and shape analyser, pycnometer
- Production and packing lab with four 20L reactors


- No handling of radioactivity => R&D cooperation

- Resin and method development “cold” => R&D partner



Organic extractant impregnated onto inert support

- « Supported Solvent Extraction » / « Solvent Impregnated Resins »
 - Distribution between two non-miscible phases
 - High density of functional groups
 - Fast kinetics/small volumes => rapid separations
 - High variety of selectivities:
 - Pure extractants, synergetic mixtures, solid extractants in diluents
- Aim: selectivity for product, no selectivity for target material
 - Elution under 'soft' conditions in small volume => labeling/injection
 - Bleeding might need to be addressed (Prefilter, AIX, CEX,...)




- Radiopharmacy/Nuclear Medicine
 - Radionuclide production
 - Cooperation with cyclotrons & reactors (NL, RN producers,...)
 - Equipment provider (targetry, synthesizer,...)
 - Separation of radionuclides from irradiated targets
 - » Diagnostics: Zr-89, Cu-64, Ga-68, Ge-68, Ti-44/5, Tc-99m, Sc-43/4...
 - **ZR Resin, CU Resin, TK200 Resin, TK400, TK201, TK202,...**
 - » Therapy: alpha emitters, Lu-177, Cu-67, Sn-117m, Sc-47...
 - **TK400, TK200, TBP Resin, CU Resin, TK211/2, TK221,...**


Domains and applications

Radiopharmacy and Nuclear Medicine


- Radiopharmacy/nuclear medicine
 - Purification of generator eluates => under development
 - Decontamination of contaminated effluents => **CL Resin**,...
 - Quality control
 - Cartridge based methods
 - **DGA sheets** (functionalized TLC, Ra-223, Ga-68, Pb-212,.... => CVUT Prague)



Radiochromatogram measured immediately after separation. Low abundant radiations of ^{227}Ac were not detected.




Radiochromatogram measured one hour after separation. Decay and ingrowth of ^{210}Pb is clearly visible.




- Original scope: Hydroxamate based resin
 - Standard for Zr separation from Y targets
 - Ready to use / no activation
 - Facile Zr elution (avoid 1M oxalic acid)
- Zr-89 production via (p,n) reaction from ^{nat}Y targets
 - High Zr/Y selectivity necessary
 - Alternative e.g. TBP Resin (=> Graves et al.)
- Application for other separations: **Ti/Sc**, **Ga/Zn**, Ge/Ga

ZR Resin – HCl

ZR Resin - D_w HCl





ZR Resin - D_w HCl



- No selectivity for Y, Sc
- High Ge/Ga selectivity at elevated HCl
- High selectivity for Zr, Ti, Nb, W over wide HCl conc. range

- No selectivity for alkalines and earth alkalines
- Lanthanides not retained
- Fe retention (dip at 2 – 3M HCl)


Zr Resin – HNO₃



- High selectivity for Zr, Ti, Nb, W over wide HNO₃ concentration range
 - Loss of selectivity at 6M HNO₃
=> Resin shows colour change
- No selectivity for Y, Sc, lanthanides, earth alkalines, most transition metals,...
- High Ge/Ga selectivity at 3M HNO₃

Zr-89 separation from Y targets

Zr separation on 1 mL ZR Resin column



- Load from 2 – 6M HCl
- Rinsing described by Holland may be used
- No activation with acetonitrile
- Quantitative Zr elution in 1.5 - 2 mL ≥ 0.05 M oxalic acid
- Clean Fe removal
- Use in commercial systems
 - Taddeo, Pinctada,...

Zr-89 separation on TBP Resin



- Method published by Graves et al.
- 400mg Y foils irradiated at 14 MeV (50 μ A)
- Dissolution in 10 mL conc. HCl
- Separation on 220 mg TBP Resin
- Load from 9.6M HCl, rinse with 20 mL 9.6M HCl
- Zr elution with 1 mL 0.1M HCl
- Zr yield: $89 \pm 3\%$, Y decontamination: 1.5×10^5
- Other applications of TBP Resin:
 - Sc isotope production (=> presentation EANM'18, Polatom)
 - Sn-117m from Cd targets
 - Sn-121m and Sn-126 determination in decommissioning samples



Nuclear Medicine and Biology
Volumes 64–65, September–October 2018, Pages 1-7




Evaluation of a chloride-based ^{89}Zr isolation strategy using a tributyl phosphate (TBP)-functionalized extraction resin

Stephen A. Graves ^a, Christopher Kutryeff ^b, Kendall E. Barrett ^b, Reinier Hernandez ^c, Paul A. Ellison ^b, Steffen Happel ^d, Eduardo Aluicio-Sarduy ^b, Todd E. Barnhart ^b, Robert J. Nickles ^b, Jonathan W. Engle ^b  

[Show more](#)


<https://doi.org/10.1016/j.nucmedbio.2018.06.003>

[Get rights and content](#)




Ti-Sc Separation (Ti-44/5)

Ti/Sc separation, ZR Resin



Zr Resin Ti/Sc





- Ti retained from (high) HCl, Sc not retained
- Ti also retained in dilute acid, Sc not => Ti generator?
- Ti elution with 0.1M citric, >0.2M oxalic acid, 0.1M H₂O₂
- Publication: Malinconico et al. use of 1M oxalic acid (J Nucl Med May 1, 2018 vol. 59 no. supplement 1 664)

JNM
The Journal of Nuclear Medicine

68Ga and 45Ti production on a GE PETtrace cyclotron using the ALCEO solid target

Mario Malinconico¹, Johan Asp², Chris Lang², Francesca Boschi¹, William Tieu², Kevin Kuan², Giacomo Guidi¹ and Prab Takhar²



➤ Ti-44 production

- 4g irradiated Sc
- 5 mL Zr Resin
- Ti-44 yield >95%
- 65.2 MBq Ti-44
- $D_f(\text{Sc}): 10^5$

Fig. 3. HCl concentration dependency of K_d for $^{44}\text{Ti}/^{46}\text{Sc}$ on ZR hydroxamate resin. Fig. 5. $^{44}\text{Ti}/^{46}\text{Sc}$ elution profile using ZR hydroxamate resin with a load of 4 g of scandium.

Use of ZR Resin as support in Ti-44/Sc-44 generators







Fig. 1. Schematic concept of a forward/reverse flow radionuclide generator.

- Direct (1 mL ZR) and reverse elution (2 mL ZR)
- 65 column volumes tested up until publication
- High Sc yields, max. Ti-44 breakthrough: $4.1 \times 10^{-4}\%$
- Obtained Sc gave labelling yields > 94%
- Generator been set-up at BNL/SBU => Poster S. Houclier ISRS 2019


Ge-68 separation from GaNi/GaCo targets

- Ge separation from GaNi or GaCo
 - Cold test on sev. 100mg GaNi per mL ZR Resin
- ZR Resin: Loading from HNO_3 , HCl or H_2SO_4
 - HNO_3 & H_2SO_4 preferred \Rightarrow GeCl_4 volatile
 - High retention/purification
- Most efficient elution e.g. NaOH or NH_4OH
- Additional purification via AIX \Rightarrow use of NH_4OH .
- Elution via 0.1M oxalic acid
- Acidification to 9M HCl then
- Conversion to 0.05M HCl via TK400 or PF



Ga-68/Ga-67 separation from Zn targets

Separation Zr resin - L: 10M HCl - E: 1.5M HCl



- Loading from:
 - dilute HNO₃ (liquid targets)
 - > 6M HCl (solid targets)
- Ga separation on ZR Resin
- Elution with 1.5M HCl
 - Too high for labeling/injection
- Ga conversion step on TK200 (TOPO)
- TK200 load from 1.5M HCl, elution in 2 – 3 BV water


Presentation EANM '17: Ga-68 from liquid targets by K. Gagnon (GEHS) et al.

Poster ISRS 2019: Ga-68 from solid targets by K. Gagnon (GEHS) et al., Tieu et al. (MITRU)

Publication Riga et al. Physica Medica 2018

- Liquid target: 1.7M ⁶⁸Zn(NO₃)₂ in 0.2M HNO₃, 12MeV, 32 min, 46 μA
- 4.3 ± 0.3 GBq EOB
- Chemical yield >75%, t = 40 min
- Purity: 99.976 ± 0.002% => Ph. Eur.

Separation TK200 resin



Cyclotron production of Ga-68

- Riga et al. Physica Medica 2018 (in press)
- Liquid target: 1.7M $^{68}\text{Zn}(\text{NO}_3)_2$ in 0.2M HNO_3
- GE PETtrace at 12MeV, 32 min, 46 μA
- 4.3 ± 0.3 GBq EOB
- Separation on ZR Resin and TK200 Resin ($t \sim 40$ min)
 - Loading of ZR Resin at $<0.1\text{M}$ HNO_3 ,
 - Rinse with 9 mL 0.1M HNO_3 .
 - Ga Elution with 5 mL 2M HCl directly onto 100 mg TK200
 - Ga Elution from TK200 with water
- Chemical yield $>75\%$,
 - 2.3 ± 0.2 GBq after separation
- Purity: $99.976 \pm 0.002\%$ => Ph. Eur.
- Target material recovery 80 – 90%
- For solid targets: single cartridge method (TK400) also under evaluation

Original paper

Production of Ga-68 with a General Electric PETtrace cyclotron by liquid target

Stefano Riga^{a,*}, Gianfranco Cicoria^b, Davide Pancaldi^a, Federico Zagni^a, Sara Vichi^c, Michele Dassenno^d, Luca Mora^e, Filippo Lodi^e, Maria Pia Morigi^d, Mario Marengo^a

S. Riga et al.

Physica B







Fig. 4. Schematic diagram of the separation process (Modular Lab, Eckert & Ziegler, Berlin).



Other examples for separations on TK200


- Zn separation from Cu



- Zn/Cu separation. Elution study, ICP-MS measurement


- Zn-65 separation. Data kindly provided by Fedor Zhuravlev, DTU

- Pt separation from Ir





- Pt/Ir separation. Elution study, ICP-MS measurement

- Sc from Ca



- Elution under rather soft conditions
- Potentially in combination with TBP Resin

- Oxime based resin
- High selectivity for Cu
 - Especially with respect to Zn and Ni
 - Widely used in mass spectrometry (Cu isotope ratio measurement)



- Load from pH >2, elution in high acid (2 – 8M)
 - Very suitable for liquid targets
 - Used for (large) solid Zn targets (=> Cu-67)
 - Loading not ideal for solid Ni targets (usually high HCl)
 - Elution in high HCl not compatible with labelling/injection
 - Evaporation or conversion
 - High purity and labelling yields

Purification of ⁶⁷Cu and Recovery of its Irradiated Zn Target

A.J. DeGraffenreid^a, R. Nidzyn^a, B. Jenkins^a, D.E. Wycoff^b, T.E. Phelps^b, A. Goldberg^a, D.G. Medvedev^a, S.S. Jurisson^b, C.S. Cutler^a

^aBrookhaven National Laboratory, C-AD/MIRP—Upton, NY (USA)

^bUniversity of Missouri, Department of Chemistry—Columbia, MO (USA)

Poster
presented at
ISRS 2017

- 13.7g Zn metal dissolved to give 312 mg ZnCl₂/mL solution at pH 2
- Loading of 60,6 mL => 18.9g ZnCL2 onto 2.4g CU Resin column => 8 mL
- Rinse with 80 mL pH2 HCl
- Eluiton in 2 x 20 mL 6M HCl
- Evaporation to dryness
- Chemical yield ~100%
- Single column D_f for Zn ~10 000
 - Additional removal indicated
- Ideally further Zn and Co removal
- Original suggestion: AIX

Nuclide	EOB Activity (mCi ± 1σ)	Cu Resin Recovery (%)			
		Load w/ Quant. Transfer	pH 2 HCl Rinse	Acid #1	Acid #2
⁶⁴ Cu	4700 ± 200	ND	ND	102	ND
⁶⁵ Zn	41.0 ± 0.8	103	ND	0.04	ND
⁵⁸ Co	63 ± 1	104	0.04	0.1	0.01

- Produced 143 mCi ⁶⁷Cu
- Quantitative recovery of radiocopper
- 99.5% radionuclidic purity—single column
- ICP-OES: 132.9 μg Cu and 1.3 mg Zn
 - Anion exchange column still needed to remove trace Zn
- Specific activity ⁶⁷Cu at EOB: 1.07 mCi/μg

Cu Resin



Robust separation that could shorten the overall processing time to separate co-produced radionuclides and large quantities of Zn from radiocopper
Cation and anion exchange columns still needed to suitably purify radiocopper

Alternatives:

- TK200: used for Zn separation from Cu targets => Cu recovered in 2M HCl
- TK201: preferred option, allows for obtaining Cu in dilute acid

Cu separation /conversion on TK201

- Cu separation from solid Ni targets and conversion usually done using AIX
 - Shrinking/swelling, trace Zn removal, elution volumes
- Under beta testing: use of TK201 (amine)
- Cu conversion: Loading onto TK201 at 6M HCl, elution in dilute acid
 - Additional Ni (& Co) removal





- Two methods for Cu-64 separation from solid Ni targets under development:
 - Direct Cu purification on TK201 for small targets and “standard” purity
 - Combined TK201/CU Resin method for high purity Cu-64
 - TK201 => direct loading onto CU Resin => pH control!


• Other on-going developments

• TK211 and TK212

- 10 μ m resin beads for improved lanthanide separation
- On-bead mix of different extractants for improved selectivity





Comparison (run under identical, non-optimum conditions):




- Increased Lu/Yb selectivity by using extractant mix of suitable ratio
- Increased radiolysis stability via alcohols

- DGA based resin for facilitated LN elution




- DGA shows excellent retention of REEs but especially heavy REE as Lu difficult to elute in 0.05M HCl
- Modification of DGA (TK221) Resin allows for faster Lu elution

- Based on Polyethylene Glycol (PEG) grafted on inert support
- Tc retention from high NaOH (preferably 5M NaOH)
- Tc retention increased by Mo
- Separation from high masses of Mo
- Elution with water
- Potential uses:
 - Radiopharmacy => direct Tc-99 production by irradiation of Mo targets
 - Decommissioning => Tc determination in decommissioning samples after sample fusion




TK202 Resin – Elution curves



- Tests show Re-Tc have similar behaviour in tested conditions
- Clean separation of Re-Tc and recovery in 5BV H₂O
- Tests at Polatom with real samples (Mo targets):
 - « Addition to the ^{99m}Tc solution even up to 3g Mo/gTK202 resin does not deteriorate distribution of ^{99m}Tc between the resin and 5M NaOH in equilibrium »
 - « ^{99m}Tc recovery from the column containing 100mg TK202 resin reached 96% »

- Scintillating Resins (PSm)
- « TK ElScint » range of products
- First: « TK TcScint »
- Developed by Uni Barcelona
- Plastic scintillator beads impregnated with selective extractants
- Direct measurement of cartridges after loading on LS counter
- Environmental monitoring => Tc-99 by LSC
- Chemical yield preferably Re/ICP-MS in effluents



Some other on-going projects

- SE Resin
 - Se-72/As-72 generator, Se-79
- n.c.a. Lu-177 separation
- Ac separation
- At separation (TK400,...)
- Functionalised polymers & silicates,...
 - e.g. DO-DGA, DE-DGA, macrocycles,...
- Improvement of radiolysis stability
- Ra separation (TK100/1, CAs)
- Cs/Rb separation (TK300)
- Impregnated scintillating beads
 - TK TcScint => Uni Barcelona
- Range of PAN based Resins
 - Decontamination of (large volumes) effluents
 - Variety of inorganic compounds embedded in organic matrix
 - Radionuclides, heavy metals,...
- Extractive discs
- Rapid tests
 - Test sticks => Uni Southampton
 - DGA Sheets (2D TLC)
 - Spin coated discs
- DGT (Diffusive Gradients in Thin Films) => 'bio-availability'
- Microfluidics



Thank you for your attention!

