



Extraction Chromatographic Resins, News Triskem International

Irtech Technical Workshop, 21-22nd May 2024 - Krakow



Overview

- Extraction chromatographic resins: principles
- Sample preparation
- Sr separations
- New Resins
 - TK-ELScint: TK-SrScint and TK-TcScint
 - TK102
 - TK221
 - TK225
- Methods under development
 - Ra separation
- Conclusions



TrisKem International



- Based in Rennes (France)
- Independent company since 02/07
 - Formerly part of Eichrom Europe
 - ISO 9001 since 2007
- Main products: extraction chromatographic resins
- Staff : 22
- R&D, QC and TechSupport group:
 - 5 RadChem PhD, 3 Technicians
- R&D: Development of new resins, techniques and applications
- Products used in several domains

Radiopharmacy
and
Nuclear Medicine


Environment and
Bioassay

Geochemistry
and
Metals Separation

Decommissioning




Radionuclides/metals separation





Extraction chromatography Principles




Organic phase impregnated onto an inert support referred as:

« Supported Solvent Extraction » / « Solvent Impregnated Resins »

Combination of liquid/liquid extraction and chromatographic techniques:

- Distribution between two non-miscible phases
- Stationary phase impregnated onto inert support (Choice of inert support depending on application (Radiolysis stability, plastic scintillators,...
- High density of functional groups
- Fast kinetics/small volumes => rapid separations
- High variety of selectivities:
 - Depending on extractant pure or in mixture used
 - **Aim: selectivity for analyte(s), no selectivity for matrix/impurities**
- Combining several cartridges can improve/facilitate separation







Types of extractants

Acidic

e.g. DIPEX (AC Resin)





Neutral

e.g. CMPO/TBP (TRU Resin), DPPP (UTEVA Resin), TK221, 'TK200'



Basic

e.g. Aliquat 336 (TEVA Resin)





Horwitz et al.

Choice of resin depends on

- Radionuclide/Metal to separate
- Matrix
- Type of measurement


=> Sample preparation = very important step for a good separation⁶

Metal Anion Complex Formation






Classic Resins: SR Resin





Sr uptake 3-8M HNO_3

Elution with low acidic HNO_3 conc or water

Ca, Na, K, stable Sr can interfere





Dietz et al. 2004





Classic Resins: DGA,N Resin

- Extractant: **DGA, Normal** (N,N,N',N'-tetra-n-octyldiglycolamide)
- Separation of Actinides, especially Am
- Robust against Fe interference, often used for soil samples





Sample preparation

Eliminate the matrix / preconcentrate the elements of interest

Non exhaustive list:

- Alkaline Fusion
- Mineralisation
- Lixiviation
- Coprecipitation
- ...



Mineralisation + Coprecipitation

Rapid Determination of Sr in 50g Soil Samples (AN-1403-10):

- Drying @ 110°C, blend/size to easily mineralise (1-2 hours)
- Add carrier* + aqua regia to mineralize the sample
- Heat to dryness + mineralise with 70% HNO₃
- Centrifuge + collect supernate (X2)
- Evaporate supernate + dissolve in 1M HCl
- Coprecipitate Sr with Ca₃(PO₄)₂
- Centrifuge + dissolve precipitate

Figure 1. Sample Preparation

Dry soil at 110°C. Blend and Size.
Remove 50g aliquot into 600mL glass beaker .

Muffle at 550°C for 1-2 hours.

Add 6mg Sr Carrier*, 50mL 70% HNO₃,
and 25mL 37% HCl.

*may need to adjust Sr carrier amount to
account for native Sr content in soil.

Heat to dryness on hot plate, medium setting.

Add 50mL 70% HNO₃. Warm sample. Transfer
solids and liquid to 250mL centrifuge tube.

Centrifuge 3500 rpm, 10 min.

Transfer supernate to 600mL beaker.

Add 25mL 70% HNO₃ to Solids. Mix and
Centrifuge. Transfer supernate to same
600mL beaker. Repeat once.

Discard solids to waste.

Evaporate supernate in 600mL beaker to dryness.

Dissolve residue in 15-20mL 1M HCl.

Transfer to 250mL centrifuge tube.

Dilute to 160mL. Add 1mL 1.25M Ca(NO₃)₂, 2mL
3.2M (NH₄)₂HPO₄, and 25mL 57% NH₄OH. Mix.

Centrifuge. Decant supernate to waste.

Continue to load solution preparation.



Mineralisation + Coprecipitation

Rapid Determination of Sr in 50g Soil Samples (AN-1403-10):

- Drying @ 110°C, blend/size to easily mineralise (1-2 hours)
 - Add carrier* + aqua regia to wet mineralize the sample
 - Heat to dryness + mineralise with 70% HNO₃
 - Centrifuge + collect supernate (x 2)
 - Evaporate supernate + dissolve in 1M HCl
 - Coprecipitate Sr with Ca₃(PO₄)₂
 - Centrifuge + dissolve precipitate
 - Coprecipitate Sr with CaF₂
 - Centrifuge + dissolve precipitate
- **Load solution**

Figure 1. Sample Preparation

Dry soil at 110°C. Blend and Size.
Remove 50g aliquot into 600mL glass beaker .

Muffle at 550°C for 1-2 hours.

Add 6mg Sr Carrier*, 50mL 70% HNO₃,
and 25mL 37% HCl.

*may need to adjust Sr carrier amount to
account for native Sr content in soil.

Heat to dryness on hot plate, medium setting.

Add 50mL 70% HNO₃. Warm sample. Transfer
solids and liquid to 250mL centrifuge tube.

Centrifuge 3500 rpm, 10 min.

Dissolve residue in 40mL 1.5M HCl.

Dilute to 170mL with H₂O.

Add 25mL 49% HF.

Mix well. Centrifuge 10 min.

Discard Supernate. ∩

Dissolve residue in 7mL 70% HNO₃,
7mL 3M HNO₃-0.25M Boric Acid,
7mL 2M Al(NO₃)₃.

Continue to load solution preparation.



Alkaline fusion + Coprecipitation

Rapid determination of SR in vegetation samples (AN-1405-10):

- Sample drying (2-4 hours)
- Wet mineralization with 70% HNO₃ + H₂O₂
- Fusion (15g NaOH for 5-10g sample) – (10 minutes)
- Dissolution + add carrier + Fe
- Coprecipitation with Ca₃(PO₄)₂ and Fe(OH)₃
- Centrifuge + dissolve precipitate
- **Load solution**

Figure 1. Sample Preparation

5-10g Vegetation sample in zirconium crucible

Muffle at 600°C.

2 hours for 5g sample.

4 hours for 10g sample.

Wet ash on hotplate with 5mL
70% HNO₃ and 5mL 30% H₂O₂.

Fuse samples with 15g NaOH at
600°C for 10 minutes.

Dissolve fusion cake with H₂O.

Transfer to 250mL centrifuge tube.

Add 125mg Fe and 4mg Sr. Dilute to 180mL.

Add 4mL 1.25M Ca(NO₃)₂, 5mL 3.2M (NH₄)₂HPO₄.

Mix. Cool in ice bath for 10min.

Centrifuge at 3500rpm. Decant Supernate.

Dissolve precipitate in 5mL warm
3M HNO₃, 7mL 70% HNO₃, and
7mL 2M Al(NO₃)₃.




Coprecipitations

Rapid Determination of ^{90}Sr in up to 40 Liter Seawater Samples (AN-1414-10):

- Acidification to pH 2 + Y+La+Fe carriers
- $\text{Fe}(\text{OH})_3$ coprecipitation
- Centrifuge
- Dissolve precipitate in 1.5M HCl + Ca carrier + HF
- CaF_2 coprecipitation
- Centrifuge + Dissolve precipitate
- **Load solution**

Figure 1. Sample Preparation






SR separations

SR separation in vegetation and soils (SR Resin)

AN-1403-10/AN-1407-10

1. Precondition SR Resin with 10ml 8M HNO_3
2. Load sample solution
3. Rinse beaker with 5ml 8M HNO_3 , add rinse to SR resin
4. Rinse SR Resin with:
 - 15ml 8M HNO_3
 - 10ml 3M HNO_3 -0.05M $\text{H}_2\text{C}_2\text{O}_4$
 - 10ml 8M HNO_3
5. Elute Sr with 20ml 0,05M HNO_3





SR separations

SR separation in vegetation and soils (SR Resin)

AN-1403-10/AN-1407-10

1. Precondition SR Resin with 10ml 8M HNO₃
2. Load sample solution
3. Rinse beaker with 5ml 8M HNO₃, add rinse to SR resin
4. Rinse SR Resin with:
 - 15ml 8M HNO₃
 - 10ml 3M HNO₃-0.05M H₂C₂O₄
 - 10ml 8M HNO₃
5. Elute Sr with 20ml 0,05M HNO₃

SR separation in 40L seawater (DGA Resin)

AN-1414-10

1. Precondition DGA,N Resin with 5ml 8M HNO₃
2. Load sample solution
3. Rinse beaker with 5ml 8M HNO₃, add rinse to DGA,N Resin
4. Rinse DGA,N Resin with:
 - 15ml 8M HNO₃ (Ca, Sr, Pb)
 - 20ml 0,05M HNO₃ (La, Ce, Sr, U)
 - 10ml 3M HNO₃-0.25M HF (U,Th)
 - 10ml 3M HCl (Ca, La, Pb)
5. Elute Y with 20ml 0,25M HCl



New Resins

TK-EIScint Resins – scintillating resins



UNIVERSITAT^{DE}
BARCELONA

=> presentation by Ines Llopart right after this presentation



TK102 Resin

- Modified version of SR Resin
 - Same crown-ether
 - Solvent, inert support and ratios => different
 - Solvent is a fluorinated alcohol
- Distribution coefficient K_d ~50% higher (Pb, Sr, Ba)
- Higher capacity (Pb and Sr)
- SR resin separation procedures can be transposed on TK102
- Specific separating methods under development



TK102 Resin – K_d values




Fig. 1: Distribution coefficients of selected elements on TK102 Resin in HNO₃
 ► Sr, Ba, Pb and Tl show high D_w in HNO₃




Fig. 2: Distribution coefficients of selected elements on TK102 Resin in HCl
 ► Pb, Tl, Sn, Sb, Ga show high D_w in HCl





Fig. 3: Distribution coefficients of Sr on TK102 Resin in 3 M HNO₃ in the presence of different salts
 ► D_w Sr decreases with KNO₃ starting at 0,05 M,
 ► no effect of NaNO₃ and Ca(NO₃)₂ up to 1 M.




TK102 Resin – Sr separation




Sr elution study in
8M HNO₃ load medium

Resins TK102 and SR
similar for the separation
of elements
Th/U/Pb/Sr/Ca/Bi/Y/Ca/Ba






TK102 Resin – Lead separation



Pb elution study with
2M HCl loading medium

Resins TK102 and SR
similar for the separation
of elements
Th/U/Pb/Sr/Ca/Bi/Y/Ca/Ba







TK221 Resin

(Papp, I., Vajda, N. & Happel, S.. *J Radioanal Nucl Chem* **331**, 3835–3846 (2022).
<https://doi.org/10.1007/s10967-022-08389-9>)

Resin based on a mixture of diglycolamide and phosphine oxide + traces long chained alcohol on inert support

Main applications in radpharm

Applications for the separation of actinides





TK221 Resin

(Papp, I., Vajda, N. & Happel, S.. *J Radioanal Nucl Chem* **331**, 3835–3846 (2022).

<https://doi.org/10.1007/s10967-022-08389-9>







Table 3 Recovery of actinide tracers from spiked water samples

	Actinides determination	
	Without Np separation	With Np separation
	Yield	Yield
	%	%
TAP water		
²³⁰ Th	90 ± 8	86 ± 7
²³⁹ Pu	108 ± 7	95 ± 7
²³⁷ Np	–	91 ± 9
²⁴¹ Am	103 ± 7	97 ± 6
²³³ U	103 ± 7	70 ± 7
SEA water		
²³⁰ Th	71 ± 7	61 ± 6
²³⁹ Pu	91 ± 7	87 ± 6
²³⁷ Np	–	93 ± 8
²⁴¹ Am	89 ± 7	92 ± 6
²³³ U	88 ± 7	59 ± 6



TK225 Resin

- TO-DGA plus ionic liquid
- Very high retention of lanthanides at medium to high acid
- Especially heavy lanthanides also very well retained at low acid concentrations
- **Main application: Removal of radiolanthanides from effluents**






Under development

Ra purification/recycling

- In case Ra needs to be purified on-column (e.g. dissolved Ra needles) => Use of TK101 for Ra retention / purification
- TK101 => crown-ether + but ionic liquid
Use of ionic liquid (=> TK101 Resin) allows for retention of Pb, Sr, Ba, Ra,... from pH ~2 – 7 without extensive extraction of other elements
- Ra separation from TK101 + further purification of Ra from Ba on TK102



Under development Ra purification/recycling




- Good Ra separation when loading from dilute HNO₃/HCl
- When eluting Ra in 3M HNO₃, Ba, Pb, Sr remain retained
- No retention of U, Th, Pt, Ir,...
- Ra eluted in 3M HNO₃
- Tl and Ba eluted in 8M HNO₃



Under development Ra purification/recycling


Further Ra fraction
purification from Ba on
TK102

Separation on 1 mL TK102 Resin (100 - 200 μ m) - ~0.5BV/min



Elution study - Ra separation from Ba on TK102 Resin in 3M HNO₃ - Ra data courtesy of N. Vajda (RadAnal)

Separation on 1 mL SR Resin (100 - 150 μ m) - 0.5BV/min



Elution study - Ra separation from Ba on SR Resin in 3M HNO₃ - Ra data courtesy of N. Vajda (RadAnal)



Conclusions

To perform the measurement of a RN/Metal from a given matrix

- Choice of resins
- Pretreatment of the sample to
 - Eliminate matrix
 - Concentrate RN/metal of interest
- Perform separation to collect RN/Metal of interest in pure fraction
- Perform measurements

Thank you for your attention !

<http://www.triskem.com>
support@triskem.fr

