



Application of SR and PB Resin in the  
determination of Pb-210 and Po-210

S.Happel – TrisKem

EVT2204728 “2023 Training Workshop of IAEA’s ALMERA  
Network on Advanced Topics in Radiochemistry Techniques:  
Lead-210 and Polonium-210”

11/10/2023



# TrisKem International



- Based in Rennes (France)
- Independent company since 02/07
  - Formerly part of Eichrom Europe
  - ISO 9001 since 2007
  - Two suppliers of SR, PB Resin (TrisKem and Eichrom) depending country
- Main product: extraction chromatographic resins
- Staff : 20
- R&D, QC and TechSupport group:
  - 4 RadChem PhD, 3 Technicians
- R&D: Development of new resins, techniques and applications
- Products used in several domains

Radiopharmacy  
and  
Nuclear Medicine

Environment and  
Bioassay


Geochemistry  
and  
Metals Separation

Decommissioning

# Extraction chromatography

## Combination of liquid-liquid extraction (LLX) and chromatography

- Liquid-liquid extraction: distribution of an element between 2 non-miscible phases
  - Most general case: organic and aqueous phase
- Chromatography: distribution of an element between a solid and a liquid phase





Change the beaker once the yellow starts to drop through.

# Extraction chromatography

## Organic extractant impregnated onto inert support

➤ « Supported Solvent Extraction » / « Solvent Impregnated Resins »

- Stationary phase impregnated onto inert support
- Distribution between two non-miscible phases
- High density of functional groups
- Fast kinetics/small volumes => rapid separations
- High variety of selectivities
- Combining several cartridges can improve/facilitate separation
- Bleeding might need to be addressed (e.g. Pb elution with 6M HCl from SR)



# Inert Support

- Choice of inert support depending on application
  - Mainly polymers
  - Radiolysis stability, plastic scintillators,...
- Inert regarding chemical reactions with
  - Organic phase/extractant
  - Elements to be isolated
  - Mobile phase
- Stable regarding mechanical actions
- High specific surface
- Defined pore and particle size distribution
- Spherical particles





# Stationary Phase

- Generally organic compounds
  - Few exceptions e.g. ABEC systems
- Non-miscible with water
- Non-volatile but soluble in volatile solvents
- Ideally fast kinetics
- High selectivity
  - Aim: selectivity for analyte, no selectivity for matrix/interferences
  - Pure extractants, synergetic mixtures, solid extractants dissolved in diluents, solids (e.g. DMG)
  - **Influence of Diluents**
- Chemically, physically and mechanically stable
- High capacity

## Types of Extractants


### Acidic

e.g. HDEHP (LN Resin)





### Neutral

e.g. CMPO/TBP (TRU Resin), DPPP (UTEVA Resin), SR & Pb Resins




### Basic

e.g. Aliquat 336 (TEVA Resin)



Metal Anion Complex Formation





Metal + Anion  $\rightleftharpoons$  Complex

Complex + Organic  $\rightleftharpoons$  Extracted



# Typically employed extractants

**Trioctylphosphine oxide  
(TOPO)**





**Octyl (phenyl)-N,N-diisobutyl-  
carbamoylmethylphosphine oxide  
(CMPO)**

**Dipentyl pentyl phosphonate  
(DPPP)**




**Bis(2-ethylhexyl) hydrogen  
phosphate  
(HDEHP)**


**Tributylphosphate  
(TBP)**



**2-Ethylhexylphosphonic acid  
mono-2-ethylhexyl ester  
(HEH[EHP])**



**N,N,N',N'-tetra-n-octyl-diglycolamide  
(TODGA)**



**4,4'(5')-di-t-butylcyclohexano-18-  
crown-6**


**Aliquat 336 (R = 8, 10)**



**Phosphine sulfide**




**Tertiary amine**




**Hydroxamate**





**Diluent:** 1-octanol




Dietz et al. 2004

- 1.0M 4,4'(5')-di-t-butylcyclohexano 18-crown-6 in 1-octanol.
- 40% (w/w) loaded onto inert chromatographic support.
- Sr max. capacity: 21 mg / 2mL resin bed,
- Advised working capacity: max. 8 mg / 2mL resin bed (Opt.: 5 mg)
- Capacity for Pb higher

# SR Resin – HNO<sub>3</sub> data

Figures 2 and 3


Acid dependency of  $k'$  for various ions at 23-25°C.  
Sr Resin



Horwitz, et al., (HP292)

Figures 4 and 5

Acid dependency of  $k'$  for various ions at 23-25°C.  
Sr Resin



Horwitz (HP199)

- Called SR Resin but higher retention of Pb
- Pb well retained at low HNO<sub>3</sub>, partial elution in water
- Attention: Po data in HNO<sub>3</sub>

# SR Resin – HCl data




Fig. 1. Distribution ratios  $K_d$  of lead, bismuth and polonium in HCl solutions.







Fig. 3 Retention of  $Po^{4+}$  on Sr resin in HCl

- High Pb and Po uptake for 2M HCl
- Pb elution with high HCl (alternative: citrate)
- To be noted: A. Knight et al showed that octanol based resins (e.g. TK400) have very similar selectivity for Po in HCl => diluent strongly involved in Po extraction

# SR Resin





- **Major interferences:** Ca, stable Sr, K, Na => lower effect on Pb
- **Monovalent interferent elimination possible:** Coprecipitation with oxalate, phosphate or iron hydroxide





- **Extractant:** 0.75M crown-ether (same crown ether as in SR Resin)
- **Diluent:** Isodecanol  
(longer carbonated chain for easier removal of Pb)
- **Pb capacity:** 29 mg.g<sup>-1</sup> PB resin (lower than SR Resin)
- **Main application:** Pb-210 in water
- Less bleeding than SR Resin (lower solubility of Isodecanol in aqueous solutions)

# Comparison PB vs SR

- Dw values PB Resin:



- Dw values SR Resin:




# Differences SR and PB Resin

- PB Resin:
  - Lower concentration (and thus amount) of crown-ether
  - Use of isodecanol instead of octanol
    - Horwitz et al showed that solutions of crown-ether in alcohols with longer chains show lower  $D_w$
    - Less bleeding
  - Deliberate choice:
    - Pb retention on PB Resin remains high
    - PB Resin allows Pb elution in water
      - Especially interesting in case of evaporation and GPC measurement
- SR Resin: typically higher Po yields
- If only Pb to be analysed PB Resin might be a good option (pricing, easier Pb elution)
- For Pb and Po (and high matrix samples) SR better choice


# New: TK102 Resin

- Modified version of SR Resin
  - Same crown-ether but diluent, inert support and ratios different
    - Higher amount of crownether
    - Short chained, fluorinated alcohol
  - Generally higher Dw and capacity for Sr, Pb, Ba
  - Improved Ra/Ba separation
  - Use of more hydrophobic diluent instead of Octanol
    - Significantly less Bleeding
- Separation methods under optimisation

Separation on 1 mL SR Resin





Separation on 1 mL TK102 Resin





# TK100/1 Resins



- Based on same crownether as SR Resin
  - Different diluents: TK100 => HDEHP, TK101 => short chained ionic liquid
  - Sr and Pb uptake also between pH 2 and 7
  - Concentration and purification on same column




- Typical applications :
  - Pb-210 in water samples (up to > 5L per 2 mL column/cartridge)
  - Sr-89/90 in water samples (up to 0,5L per 2 mL column/cartridge)
  - Sr-90 by ICP-MS (very high Zr-90 decontamination) => NPL
  - Ra-226 by ICP-MS => NPL (Load and purification in one step)
    - Agilent application note

# TK100 Resin


- Spiked water samples at pH 2
- 1L sample, 2 mL cartridge
- Load at 5 – 10 mL/min
- K and Ca not retained
  - Sr breakthrough starts at 600 mL
- Rinse:
  - 5 mL deion. water
  - 15 mL 8M HNO<sub>3</sub>
- Sr Elution with 0.5 - 2M HCl
  - Th, Pb, U remain retained
  - May be eluted sequentially
- May be used for one-column sequential separation (Sr, U, Th, Pb, Ra,...)
- Using HDEHP as diluent allows for retention of Sr, Ba, Ra, Pb from acidified water samples (e.g. pH2)
- Stability issues with Sr retention





- Use of TK101 instead => HDEHP replaced by short chained ionic liquid
- High Ra & Pb retention from dilute acid, no Th, U retention



Dw TK101 in HNO<sub>3</sub>,  
selected elements



Dw TK101 in HCl,  
selected elements




All data Russell et al NPL


- Very strong Pb retention at all HNO<sub>3</sub> concentrations and HCl < 6M HCl
- Loading of >1L on 2 mL TK101 cartridge
- Easy purification of Pb from other éléments
- Elution in high HCl (or citrate)

# TK101 Resin

### Group 1 and 2




### Group 1 and 2




All data Russell et al NPL

### Transition Metals




### Transition Metals



# Elution study TK101

- Initial testing
- 2 mL TK101 cartridge
- 5 – 10mL/min
- 1L sample (spiked tap water at pH2) – load in 10 100 mL aliquots
- Rinse with 8M HNO<sub>3</sub> and 2M HCl, Pb elution in 6M HCl (optional citrate)
- Automatisatation?



## Conclusion:

- In case of crown-ether based resins diluents have strong impact on resin selectivity
- Can be used to modify/adjust resin selectivity
- SR vs PB
  - Different crown-ether concentrations, different alcohol chain lengths
  - Facilitate Pb elution from PB Resin
  - PB Resin might be advantageous in case of Pb-210 only, especially in case of measurement by GPC (elution in water)
  - SR generally preferred in case of Pb-210 and Po-210 (higher Po yields) and in case of high matrix samples
- Using 'cation exchanging' diluents allows for direct loading of Pb (Po under examination) from acidified water samples

# Example: Pb-210 and Po-210 determination in water samples

## Context:


- European drinking water directive 98/83/EC
  - WHO recommendations
  - Calculation of Total Indicative Dose (TID)
    - $^{210}\text{Pb}$ - $^{210}\text{Po}$  activities needed
- Analysis of radioactivity in natural and drinking waters
  - Pb-210 and Po-210 part of the natural radionuclides
  - Sensitive methods needed with low DL



# Pb-210 / Po-210 in water samples

- Sample volume 1.5L (pH 2)
  - Similar methods were tested for sea water samples
- Addition of 10 mg stable Pb and Po-208 or Po-209
- Pb-Po co-precipitation with  $\text{Fe}(\text{OH})_3$ 
  - Addition of 20 mg Fe(III)
  - Heating
  - Addition of 12 mL conc. ammonium hydroxide
  - Settling / filtration or centrifugation / rinsing
- Sequential separation of Po and Pb on Sr or Pb Resin
  - Only Pb: Pb resin advantageous (pricing, facile Pb elution)
  - Pb and Po: Sr resin advantageous (better Po yields)
- Po auto-deposition on Ag or Ni disc for  $\alpha$ -spec

# Pb/Po separation on Sr or Pb resin




- After further optimisation introduced as application note/method
- Load from 2M HCl to retain Pb and Po
- Reminder, when loading from 1M  $\text{HNO}_3$  only Pb retained
- Rinse with 1M and 0.1M  $\text{HNO}_3$  to elute Po => autodeposition
- Pb elution in ammonium citrate
  - for LSC measurement
- In case of GPC measurement either
  - PB Resin => elution with water
  - SR Resin => elution in 6M HCl
    - Attn. Organics!

# Pb-210 measurement

- Pb:
  - 1 mL for AAS(/ICP-MS) measurement → determination of chemical yield
  - 9 mL for LSC measurement (or GPC)
- Pb-210 Counting:
  - Direct counting of Pb-210 via LSC ( $\varepsilon \approx 60\%$ )
    - Time dependent calibration of LSC counter
    - Preferably 2 window method
  - LSC counting after Bi-210 ingrowth
    - Better  $\varepsilon$ , easier calibration, longer waiting time
  - Also possible: evaporation of Pb eluate and GPC
    - Problem: citric acid, better elution with 20 mL dest. Water or 6M HCl

# Method development – elution studies



Alpha/Beta discrimination LSC counting

- Decontamination factors (DCF):
  - overall high for U, Th and Sr,
  - good DCF for Ra in Pb fraction,
  - rather low DCF for Ra in Po fraction - not problematic since Ra is not auto-deposited.

Interfering RN	Pb Fraction		Po Fraction	
	$\alpha$ DCF	$\beta$ DCF	$\alpha$ DCF	$\beta$ DCF
<b>Ra-226</b>	182	89	14	45
<b>U-Nat</b>	> 1000	> 1000	> 1000	> 1000
<b>Th-Nat</b>	> 1000	➤ 1000	> 1000	> 1000
<b>Sr-90</b>	NA	➤ 1000	NA	> 1000

# Po-210 autodeposition

- Evaporate eluate to near dryness ( $T < 110^{\circ}\text{C}$ )
- Convert residue into HCl form
  - Repeated evaporation with 5mL conc HCl
- Dissolve residue in 0.1M HCl
- Optional: addition of 100mg ascorbic acid
- Autodeposit Po onto Ag or Ni disc
  - T solution:  $80 - 85^{\circ}\text{C}$
  - t between 90 min and 8 h
  - stirring



# Comparison of $^{210}\text{Po}$ auto-deposition yields on Ni and Ag Discs

	Recovery (%)		FWHM (keV)	
	mean	Std (%)	mean	Std (%)
Ni (N=10)	61	30	18.4	8.2
Ag (N=10)	99	8.6	19.1	15.8

- Ag discs
  - Quantitative auto-deposition
  - Good reproducibility
- Ni discs
  - Reproducibility not stable
  - Inbetween change of supplier
- Both discs show very good resolution

# Performance data

## ■ Precision

- Repeatability  $s_r$  and Reproducibility  $s_R$  for Pb-210 and Po-210  $< 10\%$  ( $k = 1$ )

	<b>Pb-210</b>	<b>Po-210</b>
<b>Repeatability <math>s_r</math> (N=6)</b>	2.2	8.9
<b>Reproducibility <math>s_R</math> (N=12)</b>	4.4	5.2

## ■ Detection limits (DL) for $V = 1.5$ L

- Pb-210 in the order of  $20 \text{ mBq.L}^{-1}$  with counting time of 240 minutes
- Po-210 in the order of  $3 \text{ mBq.L}^{-1}$  with counting time of 1000 minutes



# Accuracy

- Validation by successful participation in intercomparison (BfS drinking water), comparison with accredited method and analysis of spiked samples

Sample ref.	Pb-210					Po-210				
	Ref. A (Bq.L <sup>-1</sup> )	UcA (Bq.L <sup>-1</sup> )	Exp. A (Bq.l <sup>-1</sup> )	UcA (Bq.L <sup>-1</sup> )	t value	Ref. A (Bq.L <sup>-1</sup> )	UcA (Bq.L <sup>-1</sup> )	Exp. A (Bq.l <sup>-1</sup> )	UcA (Bq.L <sup>-1</sup> )	t value
Water 1+	9,61 <sup>E-02</sup>	2,87 <sup>E-02</sup>	7,00 <sup>E-02</sup>	3,00 <sup>E-02</sup>	0,6	1,30 <sup>E-02</sup>	2,91 <sup>E-03</sup>	1,33 <sup>E-02</sup>	8,54 <sup>E-03</sup>	0,1
Water 2+	< LD		< LD		N/A	2,03 <sup>E-03</sup>	1,15 <sup>E-03</sup>	1,60 <sup>E-03</sup>	6,00 <sup>E-04</sup>	0,3
Spiked Water 1	1,77 <sup>E-01</sup>	1,93 <sup>E-02</sup>	1,60 <sup>E-01</sup>	3,30 <sup>E-03</sup> 3	0,9	1,60 <sup>E-01</sup>	1,93 <sup>E-02</sup>	1,39 <sup>E-01</sup>	1,84 <sup>E-02</sup>	1,0
Spiked Water 2	1,69 <sup>E-01</sup>	5,15 <sup>E-03</sup>	1,63 <sup>E-01</sup>	8,09 <sup>E-03</sup>	0,7	1,69 <sup>E-01</sup>	5,15 <sup>E-03</sup>	1,55 <sup>E-01</sup>	1,94 <sup>E-02</sup>	0,7
Spiked Water 3	1,66 <sup>E-01</sup>	5,08 <sup>E-03</sup>	1,58 <sup>E-01</sup>	8,09 <sup>E-03</sup>	0,9	1,66 <sup>E-01</sup>	5,08 <sup>E-03</sup>	1,50 <sup>E-01</sup>	2,04 <sup>E-02</sup>	0,7
Spiked Water 4	1,69 <sup>E-01</sup>	5,18 <sup>E-03</sup>	1,63 <sup>E-01</sup>	8,09 <sup>E-03</sup>	0,6	1,69 <sup>E-01</sup>	5,18 <sup>E-03</sup>	1,57 <sup>E-01</sup>	1,80 <sup>E-02</sup>	0,7
Mineralwater*	1,24 <sup>E-01</sup>	1,69 <sup>E-02</sup>	1,02 <sup>E-01</sup>	3,20 <sup>E-02</sup>	0,6	1,06 <sup>E-01</sup>	1,30 <sup>E-03</sup>	1,06 <sup>E-01</sup>	8,76 <sup>E-02</sup>	0,6

$$t = \frac{|C_{A,det} - C_{A,Ref}|}{\sqrt{U_{C_{A,det}}^2 + U_{C_{A,ref}}^2}}$$

- $t \leq 1$  • No significant difference: accurate
  - $t > 1$  • Results are not accurate
- (*t* corresponds to  $E_n$  in ISO/IEC Guide 43-1)

# Pb-210 / Po-210 in environmental and biological samples



**Dell, A.N., Curtis, M.C.: The Sequential Determination of  $^{210}\text{Pb}$  and  $^{210}\text{Po}$  in Environmental Samples using an Eichrom Resin, 6th International Conference Cogema-La Hague, 1996**

**Vajda, N et al.: A Novel technique for the Simultaneous Determination of  $^{210}\text{Pb}$  and  $^{210}\text{Po}$  Using a crown Ether, J. Environ. Radioact. 37(3), 355ff, 1997**

**Horwitz, E.P. et al.: A novel Strontium-selective extraction chromatographic resin, Solv. Ext. Ion Exch., 10(2), 1992**

**M. Villa: Isolation of lead and polonium from seawater and determination of Pb-210 by LSC, Oral presentation, LSC08, Davos 2008**

# Pb and Po in milk, meat, crop and sediment samples (Dell et al.)

- Sample volume/mass: 100mL milk, 100g crop or meat, 10g sediment
- Addition of Po-208/9 and 10 mg stable Pb
  - Preferably Po-209 (half-life,  $\alpha$ -energy)
- Wet oxidation using  $\text{HNO}_3$  and  $\text{H}_2\text{O}_2$  at  $T < 110^\circ\text{C}$ 
  - low T to avoid Po losses
  - Use of HF, « aqua regia » for soil and sediments
- Separation of Po and Pb on Sr Resin (3g)
- Po auto-deposition on Ag or Ni disc
- $\alpha$ -spec

# Elution behaviour of Bi, Pb, Po on 3g SR Resin

Vajda et al.; 1997





Fig. 2. Elution chromatogram of lead, bismuth and polonium on Sr. Spec.

- 3g columns require rather large elution volumes not very suitable for citrate elution
- Typically Pb eluted with 6M HCl => requires post-treatment

# Pb and Po in milk, meat, crop and sediment samples (Dell et al.)



- Sample preparation
- Add tracers/carrier at earliest convenience
  - 30 mg Pb in case of gravimetric yield determination
  - 1 – 10mg for spectrometry
- In case Po-210 to be determined mineralisation at  $< 110^{\circ}\text{C}$ 
  - If no Po higher temperatures may be employed
- 2 Options:
  - Complete dissolution
  - Leaching

- Complete dissolution (volumes per 5g sample):
  - 3 evaporations with 10 mL conc.  $\text{HNO}_3$  and 40 mL conc. HF
  - 3 evaporations with 30 mL conc  $\text{HNO}_3$
  - Addition of 2g  $\text{H}_3\text{BO}_3$  and 3 evaporations with 30 mL conc HCl
  - Dissolution in 30 – 50 mL 2M HCl
- Leaching (volumes per 5g sample):
  - 3 evaporations with 30 mL conc.  $\text{HNO}_3$  and 5 mL 30%  $\text{H}_2\text{O}_2$
  - 3 evaporations with 30 mL conc HCl
  - Dissolution in 30 – 50 mL 2M HCl
- Dissolution: warming for 30 min, followed by filtration

- Separation chemistry:
- 3g Sr resin column preconditioned with 2M HCl
- Load from 2M HCl
- Rinse beaker with 2 x 5 mL 5M HCl
- Rinse column with 90 mL 2M HCl (Bi removal)
- Rinse column with 60 mL 6M HNO<sub>3</sub> (Po removal)
- Elute Pb with 60 mL 6M HCl



- Preparation of the sample for LSC measurement
  - Evaporation to dryness (yellow/brown residue will remain)
  - 3 evaporations with 2 mL conc  $\text{HNO}_3$
- Option a.) Gravimetric yield determination
  - Dissolve residue in 20 mL 1M  $\text{HNO}_3$
  - Add 400 mg oxalic acid
  - Adjust pH to 6 – 7
  - Filter, dry and weight
  - Redissolve in dilute acid (1 mL 6M  $\text{HNO}_3$  or 10 mL 0.1M  $\text{HNO}_3$ )
  - Add cocktail

- Option b.) Spectrometric yield determination
  - Dissolve residue in 10 mL 0.1M HNO<sub>3</sub>
  - Withdraw aliquot for yield determination
    - In case of yield determination by ICP-MS and addition of small amounts of Pb carrier it is preferable to control Pb content of original sample to prevent over-estimation of the chemical yield
  - Add cocktail
- Measurement by LSC

- Preparation of the sample for GPC measurement
  - Addition of 2 – 3 mL conc  $H_2SO_4$ ,
  - Heating for 30 min and filtration
  - Drying of filter and weighing
  - Gravimetric yield determination
  - Filter counting after Bi ingrowth

# Po-210 via DGA Resin and MP

- Sherrod Maxwell et al.
- Rapid methods for water samples
- Co-precipitation with CaPhosphate
- Po separation via DGA Resin
- Sample preparation for alpha spectrometry by microprecipitation with Bi-Phosphate






Fig. 8 Alpha spectra for Po isotopes using bismuth phosphate (100  $\mu\text{g}$  Bi)





Thank you for your attention!

[shappel@triskem.fr](mailto:shappel@triskem.fr)




 **SUBSCRIBE TO OUR NEWSLETTER**

To keep updated with our latest developments, news and agenda for a year, subscribe to the Triskem Infos [here](#)




## Elution of a Band



Detector

Stolen from  
Bill Burnett!


# Elution experiment



Stolen from  
Bill Burnett!

# Free Column Volume (FCV)


- Volume of mobile phase in the resin bed:



- For extraction chromatography typically 0.65 mL per mL column bed
- Capacity factor  $k'$  expressed in FCV for comparison of methods using different column geometry - normalisation




# Elution curves => $k'$




$$k_1 = V_1 / \text{FCV}$$

FCV = free column volume

$k' = 5 \Rightarrow V_1 \sim 6,5 \text{ mL}$   
➤ **Need min. 15 mL for elution**




$$k_2 = V_2 / \text{FCV}$$



$$k_3 = V_3 / \text{FCV}$$

**For  $k' = 1000$**   
➤  **$V_3 \sim 1,3\text{L!}$**

# Dry Weigh Distribution Ratio




Horwitz




$$D_w = \frac{A_0 - A_s}{w(g)} \bigg/ \frac{A_s}{v(\text{mL})}$$

$$k' \sim D_w^{*0.5}$$

$$k' = f([acid])$$




# TRU Resin




From Horwitz (1)


Effect of Matrix Constituents on Americium Retention  
TRU Resin 2 M HNO<sub>3</sub>



Effect of Matrix Constituents on Neptunium Retention  
TRU Resin 2 M HNO<sub>3</sub>



# Resolution R




## Improvement of resolution R

- Decrease peak width  $\sigma$ : e.g. particle diameter, flow rate
- Increase difference of  $k'$  values : modification of stationary or mobile phase, valence, auxiliary complexants

Comparison of Elution Curves for  $\text{Sr}^{2+}$  for Two Particle sizes of Sr Resin  
Eluent 3.2 M  $\text{HNO}_3$ , 23-24°C

- Smaller particles, sharper elution bands
- Peak maximum corresponds to  $k'$





# Results (Vrecek et al. )

Table 1  
Survey of results for  $^{210}\text{Pb}$  and  $^{210}\text{Po}$  in IAEA-300 using procedure (a)–(d)

Procedures used in the processing of a sample of IAEA-300 reference material	$^{210}\text{Pb}$ ( $\text{Bq kg}^{-1}$ ) <sup>a</sup>	$^{210}\text{Po}$ ( $\text{Bq kg}^{-1}$ ) <sup>a</sup>
(Baltic Sea sediment)	(Pb recovery in %)	(Po recovery in %)
Total sample dissolution (a)	$352 \pm 42$ (91)	$401 \pm 40$ (12)
Dissolution after equilibration (b)	$351 \pm 42$ (93)	$313 \pm 31$ (19)
Leaching with HF and $\text{HNO}_3$ (c)	$348 \pm 41$ (88)	$337 \pm 34$ (42)
Leaching with $\text{H}_2\text{O}_2$ and $\text{HNO}_3$ (d)	$343 \pm 26$ (84)	$347 \pm 31$ (54)
Recommended value	360	340.5
Confidence interval	(339–395)	(273.6–361.0)

<sup>a</sup> Error margins derived from three replicate measurements.

Table 2  
Comparison of  $^{210}\text{Pb}$  and  $^{210}\text{Po}$  results by radiochemical procedure (d) and gamma spectrometry

Sample	$^{210}\text{Pb}$ ( $\text{Bq kg}^{-1}$ ) <sup>a</sup> (radiochemical: procedure d)	$^{210}\text{Pb}$ ( $\text{Bq kg}^{-1}$ ) <sup>b</sup> (non-destructive: direct $\gamma$ -detection)	$^{210}\text{Po}$ ( $\text{Bq kg}^{-1}$ ) <sup>a</sup> (radiochemical: procedure d)
River sediment 1	$63 \pm 8$	$60 \pm 13$	$65 \pm 7$
River sediment 2	$69 \pm 8$	$66 \pm 11$	$68 \pm 7$
River sediment 3	$141 \pm 7$	$144 \pm 6$	$145 \pm 6$
Soil near TPT 1	$126 \pm 13$	$121 \pm 12$	$124 \pm 12$
Soil near TPT 2	$110 \pm 13$	$100 \pm 10$	$111 \pm 9$
Soil near TPT 3	$75 \pm 9$	$80 \pm 8$	$107 \pm 9$
Needles near TPT	$68 \pm 11$	—	$69 \pm 17$
Grass near TPT	$36 \pm 4$	—	$32 \pm 6$
IAEA 368	$26 \pm 4$	$21 \pm 2$	$20 \pm 2$
(sediment)	<i>23 (15%)<sup>c</sup></i>	<i>23 (15%)<sup>c</sup></i>	
IAEA 307	$56 \pm 7$	$63 \pm 3$	$60 \pm 10$
(sea plant)	<i>59 (40–91)<sup>c</sup></i>	<i>59 (40–91)<sup>c</sup></i>	
IAEA 327	$63 \pm 7$	$57 \pm 2$	$57 \pm 9$
(soil)	<i>59 (8%)<sup>c</sup></i>	<i>59 (8%)<sup>c</sup></i>	

Note:  $^{210}\text{Pb}$  and  $^{210}\text{Po}$  are considered as being in equilibrium.

<sup>a</sup> Error margins derived from three replicate measurements.

<sup>b</sup> Error margin derived from the statistical variation of a single sample analysed.

<sup>c</sup> Certified values in italics.