

# Tritium extraction from organic-rich sample matrices using the $\text{HBO}_2$

Ian Croudace, Rich Marsh, Phil  
Warwick, Dave Reading


Raddec International Ltd

# Outline

- $^3\text{H}$  extraction by bomb combustion
- Introduction to the  $\text{HBO}_2$
- Combustion Procedure
- Safety Assessment
- $\text{HBO}_2$  (Mk II) developments
- User case study



$\text{HBO}_2$  (Mk I)  
(2011-2021)



$\text{HBO}_2$  (Mk II)  
(2022 onwards)


# Analytical requirements?

Instances where extraction of  $^3\text{H}$  from organic-rich samples is beneficial:

- Environmental monitoring
  - ↳ Often organic rich + low LOD requirements
- Nuclear decommissioning
  - ↳ Orphan wastes e.g. oils, rubbers, plastics
- Fusion reactor operational support
  - ↳ Heterogeneous soft wastes


# Bomb Combustion

- Enable  $^3\text{H}$  extraction via complete oxidation in an excess oxygen environment.
- Few commercial systems available for  $^3\text{H}$  extraction, typically a Parr 1121 used  
e.g. Moghissi et al., (1974)
- Limited sample size (< 10 g)
- Manual operation
- Incomplete oxidation (quenching)





# The Raddec Hyperbaric Oxidiser (HBO<sub>2</sub>)

- Designed specifically for <sup>3</sup>H extraction
- High capacity combustion vessel
- Operates at pressures ≤ 100 bar
- Optimised for organic-rich matrices
- Large samples can be combusted (typically up to 30 g)




# System components (I)



# Combustion Procedure



Sample pelletised  
or cut to size



Sample loaded into  
disposable silica  
crucible



Sample combusted



Measurement  
by LSC



Combustion water  
recovered under  
vacuum


# Safety Systems & Verification

- Independently verified hydraulic pressure test (136 bar)
- Finite Element Analysis - pressure compatibility
- Meets Directive 2014/68/EU (III B Prod + C2)
- High-volume Swagelok PRV
- Dual vessel closure interlocks
- Dual ignition interlocks








# Flow schematic (HBO<sub>2</sub> Mk I)





# HBO<sub>2</sub> Mk II Developments

- Simplified gas handling system and flow control
- Optional Hastelloy C/276 vessel
- Upgraded software interface & tablet PC
- Integrated exhaust outlet tube heating
- Re-design of service connections (ergonomic)






# Validation (I)



Gravimetric recovery testing based on cellulose pellet combustions - 20 g replicate pellets (left); incremental masses (right). Completed on HBO<sub>2</sub> (Mk I).

# Validation (II)

OBT in Wheat Intercomparison sample



# Case Study

## Environmental $^3\text{H}$ Measurement

### Canadian Nuclear Safety Commission (CNSC), Ottawa



All data & images courtesy of Nadereh St-Amant, CNSC, Ottawa. See also Marsh et al. (2017)




# Background

- Discharges from Canada's nuclear industry regulated by CNSC
- Principal releases from D<sub>2</sub>O moderated CANDU reactors e.g. Bruce, Darlington etc.
- Also <sup>3</sup>H processing, removal and research facilities e.g. Chalk River, SRB Technologies etc.

Independent discharge and environmental monitoring – compliance and reassurance

# $^3\text{H}$ analysis

- HTO in water, HT/HTO in air and OBT in foodstuffs
  - Fruits/berries, Vegetables, Meat





# System comparison (I)

## Parr 1211

- Max sample size 10 g
- Combustion process not visible and maybe quenched
- Manual combustion water recovery
- Cloudy / coloured combustion water
- Requires purification prior to counting by LSC

## HBO<sub>2</sub>

- Max sample size 30+ g
- Combustion progress is visible
- Integrated vacuum collection of combustion water
- Direct measurement of combustion water is possible

# System comparison (II)

## Parr 1211

tSIE % difference  
typically 12-24%



## HBO<sub>2</sub>

tSIE % difference  
typically 0-1%

tSIE % difference values measured for wheat samples, relative to distilled water quench values.

# OBT Environmental data

- 1.5 to 60 Bq/kg fresh weight

All results are Bq/kg fresh weight

System	HBO <sub>2</sub>						Parr 1121	
	Fruits and Berries	Vegetables	Potatoes	Beef	Chicken	Pork	Milk	Fish
Typical LOD (10 g sample)	2 Bq/kg (0.5 - 1 for 20-30 g sample)						3 Bq/kg	
Tritium Processing Facility	11.8 - 17.5	1.5	1.5	N/A	N/A	N/A	1.5	N/A
NPP1	1.5 - 2	1.5 - 2	2	2 - 56.5	1.5	2 - 21	2	1.8 - 2
NPP2	1.5	1.5	1.6 - 2.1	N/A	1.5 - 10.5	N/A	1.5	N/A
NPP3	1.5 - 2.8	1.5 - 1.6	N/A	N/A	N/A	N/A	1.5	1.5 - 15.9

All data courtesy of Nadereh St-Amant, CNSC, (Marsh et al. 2017).

# HBO<sub>2</sub> – Pyrolyser validation

Sample	HBO <sub>2</sub>	Pyrolyser	HBO <sub>2</sub> memory (%)
Cardboard	11300 ± 3000	11000 ± 2000	0.6
Cardboard	6000 ± 1000	8000 ± 2000	0.5
“Housekeeping”	600 ± 100	3100 ± 600	0.3
“Housekeeping”	80 ± 20	2700 ± 500	1.3
PVC	3000 ± 600	2600 ± 500	N/A

All data courtesy of Natasha Gotts, CCFE (Marsh et al. 2017).

- Good agreement for Cardboard and PVC samples
- Minimal observed memory
- Discrepancy between “Housekeeping” data are associated with the highly heterogeneous nature of the samples

# Conclusions

- The HBO<sub>2</sub> enables rapid and efficient recovery of <sup>3</sup>H from organic rich samples up to 30 g
- The system is applicable to a wide range of sample types including biota and soft-waste
- Also applicable to orphan wastes such as oils and sludge's
- Can offer improved LODs compared to thermal oxidisers.

# Acknowledgments

- Thank you to CNSC Ottawa and CCFE for sharing their data and experiences with the HBO<sub>2</sub> and allowing us to present some of these in this talk.