

The background of the slide features a close-up, shallow depth-of-field photograph of several laboratory pipettes and vials. The pipettes are arranged diagonally from the top left towards the bottom right. Some have colored caps (blue, red, green, pink). The vials are white with some having colored caps. The lighting is bright and warm, creating a yellowish glow. A solid red horizontal band is overlaid across the middle of the image, containing the main title and event information in white text.

# Overview and new Developments Radioanalytics

RadWorkshop 2024

Steffen Happel

09/09/2024



# Overview

---

- [TK200](#)
- [TK221](#)
- [Sr-90 in sea water](#)
- [TK100/1 \(Radium\)](#)
- [TK102](#)
- [Tc separation](#)
  - [TK201](#)
  - [TK202](#)
  - [TK200](#)
  - [TK-TcScint](#)
- [Upcoming: TK-SrScint](#)
- [TK400](#)
- [Calixarene based resins](#)
- [New impregnated membrane filters](#)
- [On-going projects](#)



# TK200 Resin

Based on TOPO extractant

High retention of actinides


- Higher than UTEVA, TEVA

Applications:

- Use for very efficient U removal from Pu (Wang et al.)
- Determination of Tc-99 in water samples (Ni Yuan et al.)
- Ga-68 production (in combination with ZR Resin)
- Actinide separation from water samples

Extracts actinides even at pH 1 - 2 (nitric acid)



- Preconcentration and purification of selected actinides on same column
- 'In the field'?





# Actinides on TK200

(all data Van Es et al.)




vanEs et al. [HNO<sub>3</sub>] M

VanEs et al. [HCl] M



# U/Th separation on TK200



Load: 3M HNO<sub>3</sub> or  $\geq 1$ L pH2 (HNO<sub>3</sub>)

Very clean U/Th separation

Oxalate instead of carbonate



# TK200 Resin - U/Pu separation

Recent publication by Huang et al.

Better U removal:  $D_f(U) > 10^9$

Additional U removal via He+NH<sub>3</sub>

Overall  $D_f(U) > 10^{13}$

Pu isotopes incl. Pu-238 via ICP-MS/MS






Figure S1. Analytical procedure for determination of plutonium isotopes (<sup>238</sup>Pu, <sup>239</sup>Pu, <sup>240</sup>Pu, <sup>241</sup>Pu)

Zhao Huang, Xiaolin Hou, Xue Zhao, Rapid and Simultaneous Determination of <sup>238</sup>Pu, <sup>239</sup>Pu, <sup>240</sup>Pu, and <sup>241</sup>Pu in Samples with High-Level Uranium Using ICP-MS/MS and Extraction Chromatography, Anal. Chem. 2023, 95, 34, 12931–12939, <https://doi.org/10.1021/acs.analchem.3c02526>




# Actinides on TK200 – Application


Elution of Th from 1g TK200 resin cartridge



Elution of Pu from 1g TK200 resin cartridge




Elution of U from 1g TK200 resin cartridge



- Data by N. Vajda (RadAnal)
- Method development on-going
- Direct load of U, Th and Pu from acidified water samples (here 900 mL on 2 mL TK200)
- Sequential separation on TK200
- Automatisations & ‘in the field’ preconcentration’?



# Actinides on TK200 – Application



- 'Automatisation & 'in the field' preconcentration'?





# TK200 – direct Pu load/separation

## One TK200 cartridge

- Preconcentration and purification

## Automized separation

Acidified water samples (1 L)

Flow rate 15 mL/min

U DF:  $10^4 - 10^5$

## LoD:



- 0.32  $\mu\text{Bq/L}$  Pu-239
- 2.00  $\mu\text{Bq/L}$  Pu-240




Talanta  
Volume 262, 1 September 2023, 124710




A novel strategy for Pu determination in water samples by automated separation in combination with direct ICP-MS/MS measurement

Yuyi Ni<sup>a</sup>  , Wenting Bu<sup>a</sup>, Ke Xiong<sup>a</sup>, Sheng Hu<sup>a</sup>, Chuting Yang<sup>a</sup>, Liguao Cao<sup>b</sup>




O<sub>2</sub>-He reaction/collision gas



- Reduced reagent usage
- Minimized labor intensity
- LODs for Pu down to  $\mu\text{Bq/L}$



# TK221 Resin





- DGA Resin:
- High Am retention, rather low U retention
- TRU Resin:
- Low Am retention, very high U retention
- For high matrix samples often both needed



# TK221 Resin

Resin based on a mixture of diglycolamide and phosphine oxide + traces long chained alcohol on inert support.

- Main applications in RadPharm:
  - Lu, Tb concentration from high acid and elution in small volume of dilute HCl
  - Ac-225 purification
- Applications for the separation of actinides





# TK221 Resin – actinide separation


## Modification of DGA Resin

- Contains TO-DGA / phosphine-oxide
- higher load and more radiolysis stable inert support
- Higher U retention than DGA
- Higher Am retention than TRU
  - Potential interest for Actinide separation, particularly Am

## Cooperation with Nora Vajda

## Method development for water samples


- TEVA/TK221 method
- Ideally later also soil and decommissioning samples
- Influence of Fe and Ca



Papp, I., Vajda, N. & Happel, S. *J Radioanal Nucl Chem* (2022). <https://doi.org/10.1007/s10967-022-08389-9>



# TK221 Resin



Th and Pu removed with TEVA.  
Am and U separation on TK221 Am elution before U

**Table 3** Recovery of actinide tracers from spiked water samples

	Actinides determination	
	Without Np separation	With Np separation
	Yield	Yield
	%	%
<b>TAP water</b>		
<sup>230</sup> Th	90 ± 8	86 ± 7
<sup>239</sup> Pu	108 ± 7	95 ± 7
<sup>237</sup> Np	–	91 ± 9
<sup>241</sup> Am	103 ± 7	97 ± 6
<sup>233</sup> U	103 ± 7	70 ± 7
<b>SEA water</b>		
<sup>230</sup> Th	71 ± 7	61 ± 6
<sup>239</sup> Pu	91 ± 7	87 ± 6
<sup>237</sup> Np	–	93 ± 8
<sup>241</sup> Am	89 ± 7	92 ± 6
<sup>233</sup> U	88 ± 7	59 ± 6

- Method tested on spiked tap and sea water samples
- High yields (88+ for U and Am)
- Analysis of IAEA-TEL-2021-03 WWOPT successful
- On-going: use for solid samples


Analyte	Target values		Measured values			Relative bias	MARB <sup>a</sup>	Z-score <sup>b</sup>	Test evaluation
	Mean activity concentration	Standard deviation (sd)	Activity concentration	Standard uncertainty	Relative standard uncertainty				
	Bq/kg	Bq/kg	Bq/kg	Bq/kg	%				
<sup>239</sup> Pu	5.93	2.27	5.09	0.24	4.7	14	25	0.37	Accepted
<sup>241</sup> Am	4.85	0.57	4.73	0.15	3.2	2.5	30	0.21	Accepted
<sup>244</sup> Cm	7.02	2	7.19	0.34	4.7	2.4	25	0.09	Accepted

<sup>a</sup>Maximum Acceptable Relative Bias

<sup>b</sup>Z = |X<sub>reported</sub> - X<sub>target</sub>| / s<sub>dtarget</sub>



# Tandem TK200/TK221



- Ling Zhang, Emilia Vassileva, Determination of ultra-trace level  $^{241}\text{Am}$  in marine sediment and seawater by combining TK200-TK221 tandem-column extraction chromatography and SF ICP-MS, Talanta, 271, 2024, 125724, <https://doi.org/10.1016/j.talanta.2024.125724>

## Conclusion:

For DGA separation, the decontamination factors (DFs) for Pu, U, and Th were calculated to be  $(1.8 \pm 0.4) \times 10^3$ ,  $(3.0 \pm 1.0) \times 10^4$ , and  $24 \pm 6$ , respectively.


For TK221 separation, the DFs for Pu, U, and Th were  $(1.5 \pm 0.2) \times 10^4$ ,  $(2.1 \pm 0.3) \times 10^4$ , and  $(1.2 \pm 0.1) \times 10^3$ , respectively.

It can be seen that the TK221 resin displayed remarkably better performance for the removal of Pu and Th than DGA resin, while with an excellent decontamination ability of U, close to DGA.



**=> Use for Ra-226/8 instead of DGA?**



# Ac separation => RP presentation



Load from 2 – 6M HNO<sub>3</sub>  
 Rinse with 6M HNO<sub>3</sub>  
 10M HCl => Bi removal and  
 0.05M HNO<sub>3</sub> (Fe removal)  
 Ac elution in 0.05M HCl  
 Po remains on TK221





# Sr-90 in seawater via DGA – KAERI method

LSC 2024

Raddec/Trsikem Jointworkshop, 18<sup>th</sup> Apr 2024, Porstmouth

## A simple and straightforward technique for analyzing radionuclides in seawater



18 Apr, 2024

Hyuncheol Kim ([hckim3@kaeri.re.kr](mailto:hckim3@kaeri.re.kr)), Gahyun Kim ([ghkim97@kaeri.re.kr](mailto:ghkim97@kaeri.re.kr))



한국원자력연구원  
Korea Atomic Energy Research Institute





# Sr-90 in seawater via DGA – KAERI method

## Materials and Apparatus



**AMP-PAN (or KNiFC-PAN)**

**DGA resin**

**2 mL column**



**SALT-100 (WITHTECH Ltd. ; South Korea)**

Eight peristaltic pump  
Flow rate: 10 – 100 mL min<sup>-1</sup>  
Applicable with 2 mL/ 5 mL column

[https://www.withtech.co.kr/en/busi/new\\_busiList\\_5.php](https://www.withtech.co.kr/en/busi/new_busiList_5.php)



Hidex Q-ARE

**Automated Radionuclide  
Extraction System**

Q uick  
A utomated  
R adionuclide  
E xtraction

The most advanced automated radionuclide extraction chromatography system dedicated to radionuclide separation from environmental, food and decommissioning samples.



Quick and easy-to-use unattended radionuclide extraction  
User friendly, intuitive and hassle free.

<https://www.hidex.com/>



# Sr-90 in seawater via DGA – KAERI method

## <sup>90</sup>Sr in seawater; procedure



Presentation by Hyuncheol Kim at the TrisKem / Raddec Workshop, 18.04.2024  
Available on our website  
Video also available upon request.



# Sr-90 in seawater via DGA – KAERI method

---




- Rapid automatized separation of Y-90 from sea water (60 kg)
  - Aim: determination of Sr-90
  - Faster than standard methods
  - Potential drawback: Sample needs to be adjusted to 2M HNO<sub>3</sub>  
=> large amounts of conc. HNO<sub>3</sub>
  - Ongoing: modification of DGA to allow loading from 0.2 – 0.5M HNO<sub>3</sub>  
=> significantly less HNO<sub>3</sub> conc needed for adjustment
- => TK227 Resin
- Similar approach for Am possible?
  - SALT / Q-ARE approach also for used Cs-134/7 via AMP-PAN



# TK100/1 Resins

Based on same crownether as SR Resin

- Different solvents, Sr and Pb uptake also between pH 2 and 7
- Concentration and purification on same column




Typical applications:

- Pb-210 in water samples (up to > 5L per 2 mL column/cartridge)
- Sr-90 by ICP-MS (very high Zr-90 decontamination) => NPL
- Ra-226 by ICP-MS => NPL (Load and purification in one step)
  - Agilent application note



# TK100 Resin

## Development for measurement of $^{90}\text{Sr}$ and $^{226}\text{Ra}$ by Russel and Van Es from NPL



### Rapid Analysis of Radium-226 in Water Samples by ICP-QQQ

Application Note  
Nuclear, environmental

[https://www.agilent.com/cs/library/applications/8800\\_ICP-MS\\_5991-8324EN\\_radium\\_analysis.pdf](https://www.agilent.com/cs/library/applications/8800_ICP-MS_5991-8324EN_radium_analysis.pdf)

#### Authors

Ben Russell<sup>1</sup>, Elsje May van Es<sup>1,2</sup>,  
Glenn Woods<sup>3</sup>, David Read<sup>1,2</sup>

1. National Physical Laboratory,  
Teddington, UK

2. Chemistry Department, University  
of Surrey, Guildford, Surrey, UK



Applied Radiation and Isotopes 126 (2017) 35–39

Contents lists available at ScienceDirect

Applied Radiation and Isotopes

journal homepage: [www.elsevier.com/locate/apradiso](http://www.elsevier.com/locate/apradiso)



Development of an optimised method for analysis of  $^{90}\text{Sr}$  in decommissioning wastes by triple quadrupole inductively coupled plasma mass spectrometry

B. Russell\*, M. García-Miranda, P. Ivanov

National Physical Laboratory, Hampton Road, Teddington, TW11 0LW, UK



TK100 contains HDEHP


⇒ very high Zr retention ⇒ high Zr decontamination

⇒ Sr elution in  $\geq 0.5\text{M HCl}$



# TK100 Resin

Development for measurement of  $^{90}\text{Sr}$  and  $^{226}\text{Ra}$  by



ICP-



TK100 contains HDEHP



⇒ very high Zr retention ⇒ high Zr decontamination

⇒ Sr elution in  $\geq 0.5\text{M}$  HCl




# TK101 - Radium

## TK101 Group 1 and 2



NPL

Nation's





- Ra retention from water/dilute acid up to ~0.5M HNO<sub>3</sub>/HCl
- At higher conc. selectivity closer to SR Resin/TK102 Resin

Data provided by  
Russel et al. (NPL)



# TK101 Transition Metals




Data provided by  
Russel et al. (NPL)  
24





# TK101 - Ra




Data provided by  
Russel et al. (NPL)

- No / extremely low selectivity for Th/U
- Very strong Pb retention => elution in high HCl or citrate



# Ra separation on TK101



Good Ra separation when loading from dilute Bi partially retained from 0.05M HNO<sub>3</sub>/HCl  
 HNO<sub>3</sub>/HCl

When eluting Ra in 3M HNO<sub>3</sub>, Ba, Pb, Sr remain retained

No retention of U, Th, Pt, Ir,...

Ra eluted in 3M HNO<sub>3</sub>






Further Ba removal via TK102 possible

Tl and Ba eluted in 8M HNO<sub>3</sub>



# On-going work: test of other crown ethers and diluents

Aim: two resins: a. improved version of TK101 and b. Ra Resin working at elevated acid concentration and elution in dilute acid

Compound	Name	Abbreviation
	Dicyclohexano-18-crown-6	DCH18C6
	4,4'(4,5')(5,5')-di-tert-butyl-cyclohexano-18-crown-6	DBDCH18C6
	Dibenzo-21-crown-7	DB21C7
	4,4'(4,5')(5,5')-di-tert-butyl-dibenzo-21-crown-7	BDB21C7
	Dibenzo-24-crown-8	DB24C8

**Presentation I. Dovhyi at NRC10**




# TK102 Resin

---

- Modified version of SR Resin
  - Same crown-ether
  - Dliuent, inert support and CE/D ratios => different
  - Higher Sr, Pb and Ba retention than SR Resin
  - **Less bleeding of organic materials**
    - Should lead to better shelf life
- Work by Illarion Dohvyi (Poster during ERA14), Marine Bas, Soumaya Khalfallah, Nora Vajda, Steffen Happel
- Originally optimisation for Ra/Ba separation




# TK102 Resin - Determination of $K_d$ values



Distribution coefficients of selected elements on TK102 Resin in  $\text{HNO}_3$

- Sr, Ba, Pb and Tl show high  $D_w$  in  $\text{HNO}_3$




Distribution coefficients of Sr on TK102 Resin in 3 M  $\text{HNO}_3$  in the presence of different salts

- $D_w$  Sr decreases by 30% with  $\text{NaNO}_3$  up to 1 M,
- no effect of  $\text{KNO}_3$  and  $\text{Ca}(\text{NO}_3)_2$  up to 0,05 M.




# TK102 Resin – Elution curves comparison vs SR Resin - Sr separation




TK102 Resin vs SR resin:  
Sr elution study in 8M HNO<sub>3</sub> load  
medium

Resins TK102 and SR similar for  
the separation of elements  
Th/U/Pb/Sr/Ca/Bi/Y/Ca and Ba






# TK102 Resin – Elution curves comparison vs SR Resin - Pb separation



TK102 Resin vs SR resin:  
Pb elution study with 2M HCl  
loading medium


Resins TK102 and SR similar  
for the separation of elements  
Th/U/Pb/Sr/Ca/Bi/Y/Ca and  
Ba






# Ra/Ba separation

Separation on 1 mL TK102 Resin (100 - 200 $\mu$ m) -  $\sim$ 0.5BV/min



Separation on 1 mL SR Resin (100 - 150  $\mu$ m) - 0.5BV/min



- SR Resin: high Ba breakthrough starts after 7 – 8 bed volumes
- TK102 Resin: significantly lower Ba breakthrough
- Suitable for Ba removal from Ra at 3M HNO<sub>3</sub>





# Tc-99 separation

---

Tc-99 (difficult to measure – DTM Radionuclide) – 100% beta emitter

TEVA resin allows for Tc separation but quantitative elution needs highly acidic medium

➤ New resins resins developed for loading from both acidic or alkaline media and elution in slightly alkaline or water

TK201 resin

TK202 resin

TK-TcScint




# TK201 Resin

Based on tertiary amine (weak Anion Exchanger)  
impregnated on inert support



Main applications:

- Tc separation
  - Environmental monitoring
  - Decommissioning
  - Similar selectivity to TEVA but easier to elute
    - Use of  $\text{NH}_4\text{OH}$  or 2M  $\text{HNO}_3$
- Cu separation
  - On-going development
  - Radiopharmacy





# TK201 – Actinides





- Re uptake between pH 2 – 5M HCl
- In  $\text{HNO}_3$  medium, Re fixed at pH 1-2

- High U and Pu uptake at high HCl
- Elution at low HCl

Data Russel et al. NPL



# TK201 – Dw values



Pu well retained at elevated  $HNO_3$

Tc well retained at  $c(HNO_3) < 1M$

Tc retention significantly higher in HCl

No Mo retention at  $> 0.7M HNO_3$



# TK201 Resin – Elution curve



Load of sample at pH 1-2 to retain Re & Tc

Interferences removed during load/rinse

Mo elution at 0.7M HNO<sub>3</sub>

Elution of Re @  $\geq 1M$  HNO<sub>3</sub> possible, for Tc >2M HNO<sub>3</sub>

Preferable elution options:  $\geq 0.2M$  NH<sub>4</sub>OH







# TK202 Resin

Polyethylene Glycol (PEG) grafted on inert support


Aqueous biphasic system (ABS)

Retention of chaotropic anions e.g;  $\text{TcO}_4^-$  in the presence of kosmotropic anions ( $\text{SO}_4^{2-}$ ,  $\text{CO}_3^{2-}$ ,  $\text{OH}^-$ ,  $\text{MoO}_4^{2-}$ , ...)

For samples rich in Mo: Tc yield > 90% for 6 – 8g Mo per g TK202



$D_w$  values for Tc, Re and Mo on TK202 Resin, at varying NaOH concentrations. Tc data taken from Cieszykowska et al.



$D_w$  values for Tc in 5M NaOH using 40 mg TK202 Resin, increasing amounts of Mo. Data taken from Cieszykowska et al.



# TK202 Resin


## Retention of Tc from concentrated NaOH medium (5 - 7M)

- Alkaline Fusion e.g. decommissioning samples
- Dissolution of Mo target
- Clean separation from other tested elements


Re can be used as internal standard

Elution in a small volume of water

- Eluat remains alkaline
- Load on CEX to neutralise medium + remove Na<sup>+</sup> THEN
- Load on aluminum oxide to remove last Mo traces + elution in 0.9% NaCl



Re/Tc separation from Mo on TK202 Resin




Re separation from selected elements on 2 mL TK202 Resin cartridge, load and rinse at 1 BV/min, elution at 0.25 BV/min.



# TK200 Resin

Rapid determination of  $^{99}\text{Tc}$  in water samples using  $\text{Ti}(\text{OH})_3\text{-TcO}_2$  co-precipitation and TK200 resin by liquid scintillation counting

Ni Yuan<sup>a</sup>, Quan An<sup>a</sup>, Shan Xing<sup>a,b</sup>, Xiongxin Dai<sup>a</sup>, Xiaolin Hou<sup>c,d</sup>, Yonggang Yang<sup>a</sup>, Yan Ma<sup>a</sup>







# TK200 Resin

Rapid determination of  $^{99}\text{Tc}$  in water samples using  $\text{Ti}(\text{OH})_3\text{-TcO}_2$  co-precipitation and TK200 resin by liquid scintillation counting

Ni Yuan<sup>a</sup>, Quan An<sup>a</sup>, Shan Xing<sup>a,b</sup>, Xiongxin Dai<sup>a</sup>, Xiaolin Hou<sup>c,d</sup>, Yonggang Yang<sup>a</sup>, Yan Ma<sup>a</sup>





Fig. 1. Schematic diagram of the chemical procedure for separating  $^{99}\text{Tc}$  from the water sample.

using TDCR-LSC



# TK-TcScint



UNIVERSITAT DE  
BARCELONA



Plastic scintillating beads impregnated with selective extractant

Developped by university of Barcelona

- García, Tarancón, Bagán
- « TK-ElScint » product line
  - 1<sup>st</sup> product: « TK-TcScint »
    - Quaternary ammonium + phase modifier (similar selectivity to TEVA)
    - Environment/decommissioning => Tc-99 by LSC
  - Coming soon: TK-SrScint
    - Sr and Pb
  - More products under development



# TK-TcScint



Direct measurement of the cartridge by LSC after loading and rinsing

- NO elution/evaporation/aliquoting => easy automatisaton

Chemical yield via Re/ICP-MS in eluates.



# TK-TcScint

## Use of TK-TcScint in aqueous/urine samples for Tc-99 determination (Garcia et al., TKI UGM Cambridge 2018)

### MOP:

2ml cartridge using Vacbox

1mg Re carrier


Precondition with 2ml 0.1M HCl

Load 10ml sample in 0.1M HCl

Rinse 4x2ml DI H<sub>2</sub>O

### Results

Recovery of Rhenium (by ICP-OES)	> 98.8 %
Recovery of <sup>99</sup> Tc (by LS):	> 98.8 %
<sup>99</sup> Tc Detection Efficiency (%):	89.5(0.6)
Background (cpm):	1.09
Quenching Parameter (SQP(E)):	787(7)





# Upcoming new product

- After TK-TcScint second product of impregnated Plastic Scintillation microspheres (PSm) line
- Based on SR/TK102 Resin crownether
- Samples available

## Coming soon: TK-SrScint

Plastic Scintillation microspheres (PSm) impregnated with a selective extractant  
Developed by Tarancón & Bagán at Universitat de Barcelona


- Based on SR Resin crownether and fluorinated alcohol used in TK102 Resin
- Selectivity similar to SR and TK102 Resin

Available as ready-to-use 2mL cartridges:


- Compatible with vacuum boxes
- Facile automatization

Direct measurement of cartridges :


- No elution/addition of LSC Cocktails
- Detection efficiency:
  - $t=0 > 85\%*$
  - $t=28 \text{ days} > 185\%*$
- Tested on milk\* and river water\* samples
- Sr yield  $\geq 85\%$ , deviation:  $< \pm 10\%$



TK-SrScint Resin



4,4'[5']-di-t-butylcyclohexano-18-crown-6




Sr-90 spectra on TK-SrScint at  $t=0$  (after Sr-90/Y-90 separation) and after 28 days (Y-90 ingrown), measured on a 300SL (Hidex)

Advantages:

- Less hands-on time
- Faster turn around time
- Less radioactively contaminated waste
- No mixed wastes

Literature:

- Baudat, E.; Gautier, C.; Bagán, H. et al. Optimization of a new radiochemical method based on extraction chromatographic resins and plastic scintillation for measurement of  $^{90}\text{Sr}$  in nuclear waste. *J Radioanal Nucl Chem.* <https://doi.org/10.1007/s10967-024-09396-8>. 2024
- \*I. Giménez, J. Rotger, E. Apellániz, H. Bagán, J. Tent, A. Rigol, A. Tarancón. A new method based on selective fluorescent polymers (PSresin) for the analysis of  $^{90}\text{Sr}$  in presence of  $^{210}\text{Pb}$  in environmental samples. *Applied Radiation and Isotopes*, Volume 199, 110879. <https://doi.org/10.1016/j.apradiso.2023.110879>. 2023.
- \*Marina Sáez-Muñoz, M.; Bagán, H.; Tarancón, A.; García, J.F.; Ortiz, J.; Carlos, S.; Martorell, S. Rapid methods for radiostromium determination in aerosol filters and vegetation in emergency situations using PS resin. *Journal of Radioanalytical and Nuclear Chemistry*, 322:1397-1408. <https://doi.org/10.1007/s10967-019-06779-0>. 2019.
- Marina Sáez-Muñoz, M.; Bagán, H.; Tarancón, A.; García, J.F.; Ortiz, J.; Martorell, S. Rapid method for radiostromium determination in milk in emergency situations using PS resin. *Journal of Radioanalytical and Nuclear Chemistry*, 315, 543-555. 2018
- H. Bagán, A. Tarancón, G. Rauret, J.F. García. Radiostromium separation and measurement in a single step using plastic scintillators plus selective extractants. Application to aqueous sample analysis. *Analytica Chimica Acta*, 686, 1-2, 50-56. 2011.



Application:

Sr-90 determination in environmental and decommissioning samples



# TK400 Resin

## Long chained alcohol



- First work by Knight et al. on Np/Pa separation

Retention only at high HCl concentration (>6M HCl), elution in low HCl, water,...


Main application: Pa separation

Also retains **Mo, Fe, Po, Ga, Nb**,... working on Sb

Higher Fe capacity than e.g. TRU Resin (~15mg Fe/g TK400)





Data provided by B. Russel (NPL)



46



# Fe/Mo/Nb separation



- Recovery of Fe/Nb/Mo from high HCl on TK400
- Majority of other elements removed during load and rinses (9M and 5M HCl)
- Fe/Nb and Mo eluted in dilute HCl => separation on ZR Resin
- Can also be used to remove Nb from Zr (e.g. stacked TK400/UTEVA) or Pu-241



# Calixarene based resins for Cs separation

---

AMP-PAN and KNiFC-PAN well suitable for Cs concentration from aqueous matrices but:

Cs elution difficult, leading to high matrix Cs containing solutions

- Use of  $\text{NH}_4\text{OH}$ ,  $\text{Sr}(\text{OH})_2$  followed by AIX and CEX

Use of calixarene based resins instead

Original work: TK300 Resin

- High Cs/Ba selectivity
- Load from water up to 1M  $\text{HNO}_3$
- Interference by K
- Low Cs capacity
- Home made calixarene => upscale too difficult





# Calixarene based resins

---

Aim: two resins

- Separation of Cs and Rb from neutral to weak acid and elution with strong acid
- Separation of Cs and Rb from elevated acid and elution with water or weak acid
- Use of ionic liquids or short-chained alcohols as diluents

Preparation of 13 test resins (PR) based on commercially available calixarenes


Characterisation of these resins

- $D_w$  values of different element on test resins in  $\text{HNO}_3$  and  $\text{HCl}$
- Influence of interfering ions (like potassium) on Cs separation
- Breakthrough and full capacities
- Elution tests for Rb and Cs separation


Work performed by Illarion Dovhyi, presented at last TKI/Raddec Workshop (18.04.2024) => visit our website for more data




# Tested Calixarenes



Calix[4]arene-bis(tert-octylbenzocrown-6,  $C_{72}H_{92}O_{12}$   
(BOBCalix)




1,3-alt-25,27-Bis(3,7-dimethyloctyl-1-oxy)calix[4]arene-benzocrown-6,  $C_{62}H_{82}O_8$   
(MAXCalix)




Calix[4]arene-bis[4-(2-ethylhexyl)benzo-crown-6],  
 $C_{72}H_{92}O_{12}$   
(BEBHCAlix)



# $D_W$ values of selected cations in $\text{HNO}_3$



Acid dependency of  $D_W$  for  $\text{Cs}^+$  on PR1-13 in  $\text{HNO}_3$




Acid dependency of  $D_W$  for  $\text{Rb}^+$  on PR 1-13 in  $\text{HNO}_3$

- Generally,  $D_W \text{ Cs} > D_W \text{ Rb}$ , other elements (Ag, Al, Ba, Bi, Ce, Co, Cu, Er, Eu, Ga, Hf, La, Lu, Mo, Nb, Nd, Ni, Pb, Re, Sb, Sc, Sn, Sr, Th, U, Y, Yb, Zn, Zr) not retained from  $\text{HNO}_3$
- Ionic liquid based resins: very high  $D_W$  from 0.01M to  $\geq 1\text{M}$   $\text{HNO}_3$
- Other test resins (non IL) low Cs/Rb extraction from low acid, maximum at 2 – 3M  $\text{HNO}_3$
- Generally strong decrease of retention at very high  $\text{HNO}_3$



# Elution tests with selected test resins


Separation on 1 mL PR8



Separation of Rb and Cs on PR8  
(loading in 0,01 M HNO<sub>3</sub>)

- Generally high selectivity for Cs and Rb, interferences well removed
- Cs elution required >5 BV 8M HNO<sub>3</sub>

Separation on 1 mL PR8



Separation of Rb and Cs on PR8  
(loading in 1 M HNO<sub>3</sub> solution)



# Under development: range of impregnated membrane filters

On-going work: development of  
impregnated membrane filters

First filters under beta testing:



- TK100 (DGT of Sr, Pb, Zn, LN in soil samples)
- TK201 (determination of Tc-99 in aqueous samples)
- CL for iodine,...
- 25mm and 47mm

Example: membrane filter for gross alpha  
measurement

pH 2, 10mL/min, typically 100mL samples

High retention of actinides

Glue on disc => alpha spec



Alpha spectrum, Am-241 & Pu-239, ~50mBq each



# Passive sampling

- TK100 discs
- Wagner et al.
- Passive sampling via DGT (Diffusive Gradient in Thin films) => 'bio-availability'
- Published: Sr and Pb isotope ratios in soil samples, Zn also possible
- On-going: Sr-90

analytical  
chemistry

pubs.acs.org/ac




Article

## Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils

Stefan Wagner, Jakob Santner, Johanna Irrgeher, Markus Puschenreiter, Steffen Happel, and Thomas Prohaska\*

Cite This: *Anal. Chem.* 2022, 94, 6338–6346

Read Online





# Some other on-going projects

---

- **Rapid tests**
  - Impregnated PSm resins
  - Range of 'Test sticks'
    - Suitable impregnated support
    - JCU => rapide isotope ratio analysis by MS (metallomics)
    - NPL
    - Uni Barcelona
- **Passive sampling (DGT)**
- **Separation of DTM**
  - SE Resin => Se-79
  - Zr-93, Fe, Mo, Nb,...
- **Decontamination**
  - PAN based materials (e.g. AMP-PAN)
- **Fate' of RN in the environment**
  - Separation methods
  - Mainly longer lived RN (=> therapy)
    - Ac-225/7, Lu-177(m), radioiodine,...
  - Quantification
- **In-field preconcentration**
  - Impregnated membranes
  - Cartridges
- **Microfluidics**
- **Other 'geometries' & 'Non-resin' separation materials**



# Thank you for your attention!



## SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Info here

