



Overview and new Developments RadPharm

RadWorkshop 2024

Steffen Happel

09/09/2024



Overview


- Nuclear Medicine / RadPharm
- Research interests
‘RadPharm’
- ZR Resin
 - Zr-89 from Y targets
 - Zr-89 via TBP/TK400
 - Ga-68 from Zn targets
- Cu-61/4/7
 - TK201 Resin
 - CU Resin
- Radiolanthanides
 - Tb-161 from Gd targets
- Ac-225
- Ra-226
- Tc-99m from Mo
- Quality Control – Sheets
- Other on-going R&D



Nuclear Medicine / Radiopharmacy

Use of radioactivity for imaging and treatment

- Internalisation of radionuclides and distribution in the body
- Accumulation e.g. in cancer cells
 - Imaging: PET (e.g. ^{18}F -FDG) and SPECT (e.g. Tc-99m)
 - Treatment: I-131 for Thyroid cancer
- Iodine => first theranostic system (Saul Hertz)



Source:
Jiang et al. 2019 doi: 10.3390/s19225019

Renewed interest in use for therapy

- Bayer acquires Algeta => Ra-223 (Xofigo)
- Novartis (AAA and Endocyte) => Lu-177
- Generally use of alpha or beta emitter
 - Less frequent: Auger-Meitner emitters

Choice of Radionuclide			
Nuclide	T _{1/2}	emission	mean path length
I-125	60.0d	auger	→ 10nm
At-211	7.2h	alpha	→ 65nm
Lu-177	6.7d	beta/gamma	→ 0.7mm
Cu-67	2.58d	beta/gamma	→ 0.7mm
I-131	8.04d	beta/gamma	→ 0.9mm*
Sm-153	1.95d	beta/gamma	→ 1.2mm
Re-186	3.8d	beta/gamma	→ 1.8mm
P-32	14.3d	beta	→ 2.9mm
Re-188	17h	beta/gamma	→ 3.5mm
In-114m	50d	beta/gamma	→ 3.6mm
Y-90	2.67	beta	→ 3.9mm*

*¹³¹I ----> 3mm dia. ⁹⁰Y ----> 2cm dia.


The University of Nottingham | Meldon et. al. Radiother. Oncol. 1998





Thera(g)nostics

« Treat what you see and see what you treat »

- Step towards personalized medicine
- Injection of targeted radiotracer
 - Labelling with positron or gamma emitter for imaging
 - Size / position of tumor, tracer distribution in body
 - Ideally allows dose calculation/adjustment
 - Decision whether treatment is suitable upfront
 - Labelling with alpha, beta (or Auger) emitter for treatment
- Theranostic pair for imaging and treatment
 - ‘Real’ theranostic pair: Cu-64/7, Sc-44/7, Tb, Pb,...
 - Other theranostic pairs: Ga-68/Lu-177, Ga-68/Ac-225,...
 - Sufficiently similar chemistry



Source: <https://wilson.chem.cornell.edu/research/>



Source: <https://uihc.org/health-topics/what-theranostics>



Most promising systems: PSMA & FAPI

- PSMA: Treatment of metastatic castration resistant prostate cancer
- PSMA-617, PSMA-11, ...
- Ga-68 and Lu-177 or Ac-225
- i.e. Vision study (Novartis)
- Large interest in FAPI
 - Detection of 28 different cancers






Figure 1. Structures of (a) PSMA-617 and (b) PSMA-11.

Source: Eppard et al. 2017 doi: 10.7150/thno.20586



Source: Yadav et al. 2018 doi:10.1007/s00259-016-3481-7




SNMMI image of the year 2019
Kratochwil et al. doi: 0.2967/jnumed.119.227967



Radionuclide production

- 'Legacy materials'
 - "Th"/Pb-212, Th-229/Ac-225
- Cyclotron
 - Irradiation of targets e.g. with protons (i.e. $^{68}\text{Zn}(p,n)^{68}\text{Ga}$)
- Reactors (or other neutron sources)
 - Fission (e.g. Mo-99)
 - « Neutron reactions »
 - 'Carrier added' Lu-177 \Rightarrow Lu-176 (n, γ) Lu-177
 - 'Non-carrier added' Lu-177 \Rightarrow Yb-176 (n, γ) Yb-177 \rightarrow Lu-177 + β^-
- Common challenges:
 - Large excess of matrix (target material)
 - Very high decontamination factors required
 - Especially cyclotron produced radionuclides:
 - often quite short half-life of product
 - Very high radioactivity levels \Rightarrow increased radiation stability






Research interests - Radiopharmacy

• Radionuclide production/purification

- Resin and method development 'cold'
 - Cooperation with cyclotrons & reactors (NL, RN producers,...)
 - Equipment provider (targetry, synthesizer,...)
- Separation of radionuclides from irradiated targets
 - Diagnostics: Zr-89, Cu-64, Ga-68, Ge-68, Ti-44/5, Tc-99m, Sc-43/4...
 - Therapy: alpha emitters, Lu-177, Tb-161, Cu-67, Sn-117m, Sc-47...
- Requirements for resins:
 - Choice of right resin particularly important
 - No selectivity for target material, high selectivity for product
 - Elution under 'soft' conditions in small volume => labelling/injection
 - Fast kinetics
 - Combining several resins can facilitate the separation
 - Conversion (high acid to dilute acid)
 - Removal of impurities upfront

Radiopharmacy
and
Nuclear Medicine





Research interests - Radiopharmacy

Radiopharmacy
and
Nuclear Medicine

- Quality control
 - Cartridge based methods (e.g. Sr-90 in Y-90,...)
 - Use of “TK-ELScint cartridges”?
 - **“Sheets”**
 - p.ex. DGA sheets (functionalized TLC for Ra-223, Ga-68, Pb-212,... => CVUT Prague), CU Sheets,...
- Decontamination of effluents/waste (Ge-68, lanthanides, radioiodine,...)
- ‘Recycling’/valorization of long-lived RNs (Ge-68,...) and target materials
- Radiolysis stability (polymer, radical scavengers,...)
- **Determination of radionuclides (mainly used in therapy, generally Lu-177 and Ac-225) in environmental and bioassay samples**



ZR Resin

Original scope: Hydroxamate based resin

- Different from Holland et al.
- Standard for Zr separation from Y targets
- Ready to use / no activation
- Facile Zr elution (avoid 1M oxalic acid)

Zr-89 production via (p,n) reaction from ^{nat}Y targets

- High Zr/Y selectivity necessary
- Alternative e.g. TBP Resin (=> Graves et al.)


Application for other separations: **Ti/Sc, Ga/Zn, Ge/Ga**

On-going work => improvement of radiolysis stability




ZR Resin – HCl

ZR Resin - D_w HCl




ZR Resin - D_w HCl



- No selectivity for Y, Sc
- High selectivity for Zr, Ti, Nb, W over wide HCl conc. range
- High Ge/Ga selectivity at elevated HCl
- No selectivity for alkalines and earth alkalines
- Lanthanides not retained
- Fe retention (dip at 2 – 3M HCl)



Zr Resin – HNO₃




- High selectivity for Zr, Ti, Nb, W over wide HNO₃ concentration range
 - Loss of selectivity at 6M HNO₃
=> Resin shows colour change
- No selectivity for Y, Sc, lanthanides, earth alkalines, most transition metals,...
- High Ge/Ga selectivity at 3M HNO₃



Zr-89 separation from Y targets

Zr separation on 1 mL ZR Resin column




- Load from 2 – 6M HCl
- Rinsing described by Holland may be used
- No activation with acetonitrile
- Quantitative Zr elution in 1.5 - 2 mL ≥ 0.05 M oxalic acid
- Clean Fe removal
- Use in commercial systems
 - Taddeo, Pinctada,...



Zr-89 separation on TBP Resin


- Frequent request: Zr elution without oxalate
- Method published by Graves et al.
 - 400mg Y foils irradiated at 14 MeV (50 μ A)
 - Separation on 220 mg TBP Resin
 - Load from 9.6M HCl, rinse with 20 mL 9.6M HCl
 - Zr elution with 1 mL 0.1M HCl
- Zr yield: $89 \pm 3\%$, Y decontamination: 1.5×10^5
- Zr elution should also be possible with citrate, phosphate, oxalate...
- (Fe and) Nb removal not ideal
- Improvement => use of TK400 Resin
- Implemented successfully by Lyashchenko et al.
- Two TK400 cartridges for high Fe removal



Nuclear Medicine and Biology
Volumes 64–65, September–October 2018, Pages 1–7



Evaluation of a chloride-based ^{89}Zr isolation strategy using a tributyl phosphate (TBP)-functionalized extraction resin

Stephen A. Graves ^a, Christopher Kutryreff ^b, Kendall E. Barrett ^b, Reinier Hernandez ^c, Paul A. Ellison ^b, Steffen Happel ^d, Eduardo Aluicio-Sarduy ^b, Todd E. Barnhart ^b, Robert J. Nickles ^b, Jonathan W. Engle ^b 

[Show more](#)

<https://doi.org/10.1016/j.nucmedbio.2018.06.003>

[Get rights and content](#)

Nuclear Medicine and Biology 136–137 (2024) 108943



Contents lists available at [ScienceDirect](#)

Nuclear Medicine and Biology

journal homepage: www.elsevier.com/locate/nucmedbio



^{89}Zr]ZrCl₄ for direct radiolabeling of DOTA-based precursors[☆]

Serge K. Lyashchenko ^{a,b,*}, Tuan Tran ^a, Steffen Happel ^c, Hijin Park ^a, David Bauer ^b, Kali Jones ^b, Tullio V. Esposito ^b, NagaVaraKishore Pillarsetty ^b, Jason S. Lewis ^{a,b,d}

^a Radiochemistry and Molecular Imaging Probe Core Facility, Memorial Sloan Kettering Cancer Center, New York, NY, USA

^b Department of Radiology, Memorial Sloan Kettering Cancer Center, New York, NY, USA

^c Triskem International, Inc., USA


^d Program in Molecular Pharmacology, Memorial Sloan Kettering Cancer Center, New York, NY, USA





Use of TK400 for Fe/Nb removal


Elution study TBP Resin



- On TBP only: Fe and Nb follow Zr
- Removal of Fe & Nb upfront possible using TK400 Resin
- Test with stacked 2 mL TK400/TBP cartridges

- Load and Rinse at 10M HCl with TK400 stacked above TBP
- Splitting of cartridges and separate elution with dilute Hcl

Elution study TK400/TBP Resin




- TBP => ZR only
- TK400 => Fe & Nb
- For Nb/Fe separation => Fe(II)
- Y passes through both
- Potential for Nb separation from Zr targets




Ga-68 separation from Zn targets

- Irradiation of Zn-68 targets in cyclotron
- Ga-68 separation on ZR Resin
 - No selectivity for Zn (target material)
 - Loading possible from:
 - dilute acid (**liquid targets => typically HNO₃**)
 - >6M HCl (**solid targets**)
 - Rinse under loading condition
 - Elution with ~1 - 2M HCl
- Too acidic for injection or labelling



⇒ **New IAEA TechDoc:**

<https://www-pub.iaea.org/books/IAEABooks/13484/Gallium-68-Cyclotron-Production>¹⁵



- Ga-68 'conversion' necessary
 - Evaporation & dissolution difficult to automatize
- Easier => use of another resin
- TK200 Resin (TOPO) load from 1 - 2M HCl
- Rinse with e.g. 1 - 2M HCl
- Elution in 2 – 3 BV water, dilute acid,..



Cyclotron production of Ga-68

Rodnick et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:25
<https://doi.org/10.1186/s41181-020-00106-9>

EJNMMI Radiopharmacy
and Chemistry

RESEARCH ARTICLE

Open Access

Cyclotron-based production of ^{68}Ga , $[^{68}\text{Ga}]\text{GaCl}_3$, and $[^{68}\text{Ga}]\text{Ga-PSMA-11}$ from a liquid target

Melissa E. Rodnick¹, Carina Sollert², Daniela Stark³, Mara Clark¹, Andrew Katsifis³, Brian G. Hockley¹, D. Christian Parr², Jens Frigeli², Bradford D. Henderson¹, Monica Abghari-Gerst¹, Morand R. Piert¹, Michael J. Fulham⁴, Stefan Eberl^{5*}, Katherine Gagnon^{2*} and Peter J. H. Scott^{1*}




Fig. 4 Two-column approach for ^{68}Ga chemical separation

Table 1 High level schemes of $[^{68}\text{Ga}]\text{GaCl}_3$ purifications

	Scheme A*	Scheme B
① ZR Load	< 0.1 M HNO_3	
② ZR Wash	15 mL 0.1 M HNO_3	
③ ZR Elution / Trapping on TK200	5–6 mL ~ 1.75 M HCl	
④ TK Wash	–	3.5 mL 2.0 M NaCl in 0.13 M HCl
⑤ TK Elution	H_2O	1–2 mL H_2O followed by dilute HCl to formulate

*Process as reported previously (Nair et al. 2017)

J. Kumlin et al.

ZR, LN & TK200 for solid targets

- High Ga-68 activities
- ARTMS/Odense: 10 Ci production

One column separation possible using TK400 Resin => solid targets

- Ga retention on TK400 from high HCl
- No Zn retention
- Faster kinetics than ZR Resin

W. Tieu et al. use of single TK400 cartridge for solid Zn targets

Svedjehed et al. use of TK400/A8/TK200 for solid Zn targets

Demystifying solid targets: Simple and rapid distribution-scale production of $[^{68}\text{Ga}]\text{GaCl}_3$ and $[^{68}\text{Ga}]\text{Ga-PSMA-11}$


Johan Svedjehed, Martin Pärnaste, Katherine Gagnon*

Cyclotrons and TRACERcenter, GEMS PET Systems AB, GE Healthcare, Uppsala, Sweden




Ongoing work on TK200

- Pd separation from Rh




- Pd/Rh separation. Elution study, ICP-MS measurement

- Pt separation from Ir




- Pt/Ir separation. Elution study, ICP-MS measurement

- Zn separation from Cu





- Zn/Cu separation. Elution study, ICP-MS measurement



- Zn-65 separation. Data kindly by Fedor Zhuravlev, DTU



Ru separation on TK200






Ag/Pd separation

RESEARCH ARTICLE

Open Access


Chromatographic separation of silver-111 from neutron-irradiated palladium target: toward direct labeling of radiotracers

Marianna Tosato^{1,2}, Andrea Gandini³, Steffen Happel⁴, Marine Bas⁴, Antonietta Donzella^{5,6}, Aldo Zenoni^{5,6}, Andrea Salvini³, Alberto Andrichetto⁷, Valerio Di Marco² and Mattia Asti^{1*}





Ti-Sc Separation (Ti-44/5)



- Ti retained from (high) HCl, Sc not retained
- Ti also retained in dilute acid, Sc not => Ti generator?
- Ti elution with 0.1M citric, >0.2M oxalic acid, 0.1M H₂O₂
- K. Olguin: https://www.triskem-international.com/scripts/files/5fc95b3398a614.03970900/olguin_sfu_vugm20_production-and-purification-of-titanium-45-for-positron-emission-tomography.pdf
- Publication:
 - Malinconico et al.: J Nucl Med May 1, 2018 vol. 59 (supplement 1 664)

JNM
The Journal of Nuclear Medicine


68Ga and 45Ti production on a GE PETtrace cyclotron using the ALCEO solid target

Mario Malinconico¹, Johan Asp², Chris Lang², Francesca Boschi¹, William Tieu², Kevin Kuan², Giacomo Guidi¹ and Prab Takhar²



Cu-61/4 separation on TK201

- Cu-64 separation from solid Ni-64 targets
 - Target dissolution in high HCl
 - Load and rinse at 6M HCl
 - Ni removal and recovery/recycling
 - Co elution with 4 – 5M HCl => Co separation
 - Cu elution with 0.5M HCl
 - Zn remains retained (Ga and Fe partially co-elute)
=> requires further treatment



- Improvements:
 - Preferred alternative: Use of TBP (or TK400) upfront for Fe/Ga removal
=> allows for Cu elution in 0.05M HCl

Svedjehed et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:21
<https://doi.org/10.1186/s41181-020-00108-7>

(2020) 5:21

EJNMMI Radiopharmacy and Chemistry


RESEARCH ARTICLE

Open Access

Automated, cassette-based isolation and formulation of high-purity [⁶¹Cu]CuCl₂ from solid Ni targets



Johan Svedjehed¹, Christopher J. Kutyriff², Jonathan W. Engle^{2,3} and Katherine Gagnon^{1*}



- Gagnon et al. use of NaCl/HCl for better pH control of eluate
- Also be used for Zn separation 21





CU Resin

TK201 can not be used for Cu separation from Zn targets (e.g. Cu-67)

Use of oxime based CU Resin instead

High selectivity for Cu particularly with respect to Zn, Ni, Fe, Co,...



Load from pH >2, elution in high mineral acid (2 – 8M)

- Used for (large) solid **Zn** targets (=> Cu-67)
- Not ideal for solid Ni targets (usually high HCl) => TK201
 - Works for liquid targets (pH 2 – 3) => Fonseca et al.
- Elution in high HCl not compatible with labelling/injection
 - Evaporation/redissolution or
 - Conversion to dilute HCl e.g. via TK201 (additional Zn removal) e.g. Kawabata et al.



Cu-67 at BNL (DeGraffenreid et al.)

Purification of ^{67}Cu and Recovery of its Irradiated Zn Target

A.J. DeGraffenreid^a, R. Nidzyn^a, B. Jenkins^a, D.E. Wycoff^b, T.E. Phelps^b, A. Goldberg^a, D.G. Medvedev^a, S.S. Jurisson^b, C.S. Cutler^a

^aBrookhaven National Laboratory, C-AD/MIRP—Upton, NY (USA)

^bUniversity of Missouri, Department of Chemistry—Columbia, MO (USA)

Poster presented
at ISRS 2017

Cu Resin

Recovery (%)

Nuclide	EOB Activity (mCi \pm 1 σ)	Load w/ Quant. Transfer	pH 2 HCl		
			Rinse	Acid #1	Acid #2
^{64}Cu	4700 \pm 200	ND	ND	102	ND
^{65}Zn	41.0 \pm 0.8	103	ND	0.04	ND
^{58}Co	63 \pm 1	104	0.04	0.1	0.01

- Produced 143 mCi ^{67}Cu
- Quantitative recovery of radiocopper
- 99.5% radionuclidic purity—single column
- ICP-OES: 132.9 μg Cu and 1.3 mg Zn
 - Anion exchange column still needed to remove trace Zn
- Specific activity ^{67}Cu at EOB: 1.07 mCi/ μg

Cu Resin

Robust separation that could shorten the overall processing time to separate co-produced radionuclides and large quantities of Zn from radiocopper
Cation and anion exchange columns still needed to suitably purify radiocopper

Alternatives to AIX- use of TK201:

- 13.7g Zn metal dissolved to give 312 mg ZnCl_2 /mL solution at pH 2
- Loading of 60,6 mL \Rightarrow 18.9g ZnCl_2 onto 2.4g CU Resin column \Rightarrow 8 mL
- Rinse with 80 mL pH2 HCl
- Elution in 2 x 20 mL 6M HCl
- Evaporation to dryness
- Chemical yield \sim 100%
- Single column D_f for Zn \sim 10 000
 - Additional removal indicated
- Ideally further Zn and Co removal
- Original suggestion: AIX

- Cu elution with 6M HCl directly onto TK201
- Cu elution from TK201 in dilute acid
- Optional: rinse with NaCl/HCl for better pH control




Lu-177/Tb-161

Lu-177 still more frequently used but Tb-161 getting strong interest

- Part of the 'Swiss knife of nuclear medicine' => Tb isotopes

Similar production for both



Tb 149		Tb 152		Tb 155	Tb 161
4.2m	4.1h	4.2m	17.5h	5.32d	6.90d
ε	ε	γ283;	ε	ε	β ⁻ 0.5; 0.6...
β ⁺	α3.97	160...	β ⁺ 2.8...	γ87;	γ 26; 49; 75...
α3.99	β ⁺ 1.8	ε: β ⁺ ...	γ 344;	105;...	e ⁻
γ796;	γ352;	γ344;	586;	180, 262	
165...	165...	411...	271...		

Terbium: a new 'Swiss army knife' for nuclear medicine

Source: <https://cerncourier.com/a/terbium-a-new-swiss-army-knife-for-nuclear-medicine/>

- Irradiation of several hundreds of mg or more
- Upscale on-going (incl. recycling) => typically 1g


Prepacked PP columns available

- 4cm x 30cm (375 mL), 2.5cm x 30cm, 1.5cm x 30cm & 1.1cm x 30cm
- Connection: ¼" 28G, up to ~10bar
- QC/CoA per column (peak asymmetry) for TK211/2/3
- TK221 => dry packing







Lanthanide separation on TK211/2/3



HDEHP (LN)




HEH[EHP] (LN2)



H[TMPeP] (LN3)


Cyanex 572



- Mixtures of different extractants
- Optimized for high radiation stability



Simplified method for Lu separation from 500 mg Yb – TK211/2 & TK221



Sequential separation step (direct load from TK212 onto TK211 for polish)

Unfortunately complete sequential TK213=>TK212=>TK211 didn't work out

Can be upscaled (larger columns,...)

Further optimisation on-going



Tb separation from 1000 mg Gd targets

Irradiated target typically oxide => dissolved in $>3M$ HNO_3



- For separation solution needs to be dilute acid

Conversion via TK221 Resin

Sequential separation on TK212/TK211

Final conversion to dilute HCl on TK221 + trace nitrate removal on AIX


Mainly Tb-161, also Tb-155





Tb separation from 1000 mg Gd targets




- Initial separation on TK212 – 150 mL column (30cm x 2.5cm)
- Allows for working with up to 2g of Gd
- Gd recovery => very expensive & difficult to find
- Tb separation from Gd and Dy – ideally using online detection
- Fine purification on TK211 (29 mL)



Tb separation from 1000 mg Gd on TK212 (147 mL column)



Increased amounts of Gd



On the same TK212 column

More Gd => earlier elution

- At 3g Gd start of breakthrough
- More than 3g possible? Tb needs to remain retained...

Little effect on Tb


Small impact on Dy

Tb / Dy separation remains good.

=> Further improvement via EtO addition?



Tb purification on TK211



- Direct load of Tb fraction from TK212 onto TK211 (29 mL – 30cm x 1.1cm)
- Gd breakthrough during load & rinse with 0.5M HNO_3 (alternatively HCl)
- Tb elution (Dy sufficiently well removed before) preferably in **>3M HNO_3**
- Conversion to dilute HCl via TK221, A8 for nitrate removal



TK221 Resin

DGA well suited for 'conversion' and purification (Ca, Al, Fe,... removal)

- Convert Lu from high nitric acid to dilute HCl


Elution of heavy lanthanides at as low volumes as possible

- small volume preferred => high activity concentration


TK221 Resin

- DGA / phosphine-oxide, improved radiolysis stability (inert support, scavenger,...)
- Better La and U retention
- Lu & Tb eluted in small volume in dilute HCl => drawback, no group RE separation possible

DGA Resin



TK221 Resin



- **TK222 Resin => DGA, B/phosphine oxide**





TK225 Resin

TO-DGA plus ionic liquid => originally failed experiment

Very high retention of lanthanides at medium to high acid

Especially heavy lanthanides also very well retained at low acid concentrations

Main application: Removal of radiolanthanides from effluents





Ac-225 separation

Ac-225 separation chemistry well established

- DGA (mainly B is used) allows for facile Ac/Ra and La/Ac separation
- Problem: availability of branched and normal DGA Resin
 - Eichrom/Northstar situation & Eichrom/TrisKem contract
 - Imperfect La/Ac separation
 - Radiolysis stability sufficient?

Ongoing beta testing:

- Use of TK221 (TO-DGA / phosphine oxide) or
- TK222 (TEH-DGA / phosphine oxide)
 - Focus on La/Ac separation
 - Ac elution
 - Resalting possible? Ac nitrate => Ac chloride
 - Improved radiolysis stability?



TK221/2 Resins – Ac Dw values

- Ac Dw data determination ongoing
- Work with several groups
- Upcoming publication

CONFIDENTIAL

Data courtesy of N. Vajda (RadAnal)





FIG. 1 The D_W values for ^{225}Ac between TK221 and nitric acid




FIG. 2 The D_W values for ^{225}Ac between TK221 and TK222 and hydrochloric acid


Data courtesy of
O. Lebeda (UjV
Rez)
Upcoming
publication



Ac-225 separation

Two TK221 cartridges for removal of impurities incl. La

- In case La can be excluded step 2 only




Step 1 TK221:

Target dissolved in 2 – 4M HNO₃

Ra, Ba, Pb, Sr,... removal with 4M HNO₃

Ac elution in 14M HNO₃ (LNs retained)



Step 2 TK221:

2x diluted eluate from first TK221

Rinse with 6M HNO₃ and optional rinses with:


10M HCl => Bi removal and

0.05M HNO₃ (Fe, Po removal) => can be inverted

Ac elution in 0.05M HCl




TK221 Resin – Ac separation – step one



- In case LN need to be removed
- Two step procedure

First Ac / LN separation TK221

- Load from elevated HNO₃
- Ac elution in very high HNO₃
- LNs, U, Th retained
- Particular attention to Pb/Sr
 - Elution in 4M HNO₃





TK222

- Load from elevated HNO₃
- Ac elution in very high HNO₃
- LNs, U, Th retained
- Particular attention to Pb and Sr
 - Pb Elution in 4M HNO₃
 - Sr elution in 12M HNO₃

TK221 preferred option



TK221 Resin – Ac separation – step two



- Second separation
- TK221
 - Dilute x2 => load
 - Bi removal 10M HCl
 - Fe removal in 0.05M HNO₃
 - Ac elution in 0.05M HCl
 - Important: Lanthanides need to be removed upfront (1st TK221)
 - Additional purification on TK101 possible (Ra, Ba, Pb, Sr)
- Alternative: TK222
 - Sharper Ac elution
 - No rinse with 0.05M HNO₃
 - Only in case of absence of Fe
- TK221 preferred in case of presence of Fe

Ac elution from TK221 and TK222




FIG. 3 Elution profile of ^{225}Ac loaded on the TK221 column in 3M nitric acid, rinsed with 4M nitric acid and eluted into 0.05M hydrochloric acid




FIG. 4 Elution profile of ^{225}Ac loaded on the TK222 column in 3M nitric acid, rinsed with 4M nitric acid and eluted into 0.05M hydrochloric acid




FIG. 5 Re-salting of ^{225}Ac loaded on the TK221 column in 0.05M nitric acid and eluted into 0.05M hydrochloric acid


- Sharp Ac elution from TK221, even sharper from TK222
- ‘Resalting’ from Ac nitrate to chloride form possible on TK221 (not TK222)
 - Load from 0.05M HNO_3 , elution in 0.05M HCl

CONFIDENTIAL




Ac-225 separation – Optional additional purification

Optional:
Pb removal on TK102



Optional Pb removal step (TK102)
Eluate of step 1 dilute by x2
Load through TK102
Pb and Sr retained, Ac passes through



Optional:
Ra, Pb, Sr, Ba removal on TK101



Optional Pb, Ra, Sr,... removal step
(TK101)
Pass Ac fraction (0.05M HCl) through
TK101
Ac passes through, Ra, Pb, Sr,... retained



Optional: TK101 purification step



Optional:
Ra, Pb, Sr, Ba removal on TK101

Load:

Ac-225 in <5 mL 0.05M HCl

1

2

Rinse:

1 mL 0.05M HCl



Ac-225 in 0.05M HCl
Ra, Pb, Sr remain retained

Data courtesy of B. Russel (NPL)

Optional Pb, Bi, Ra, Sr,... removal step
(TK101)

Pass Ac fraction (0.05M HCl) through TK101
Ac passes - Ra, Pb, Sr, Bi,... retained



Ra purification / recycling

Needles and other Ra sources often contain Pt, Ir, Au, Ba besides Ra.

Ra generally present as RaSO_4

Suggestion: Work-up following Matyskin et al.

Destruction of the needle, generally cutting (higher losses) or dissolution in aqua regia

Conversion of $\text{Ra}(\text{Ba})\text{SO}_4$ via heating with Na_2CO_3 solution

=> Load onto TK101 for purification



Alternatively: dissolution in EDTA solution

=> TK101 allows for Ra extraction from EDTA at pH 4




TK101 - Radium

TK101 Group 1 and 2



NPL

Nation:




- Ra retention from water/dilute acid up to $\sim 0.5\text{M HNO}_3/\text{HCl}$
- At higher conc. selectivity closer to SR Resin/TK102 Resin

Data provided by
Russel et al. (NPL)



Ra separation on TK101



Good Ra separation when loading from dilute Bi partially retained from 0.05M HNO₃/HCl
 HNO₃/HCl

When eluting Ra in 3M HNO₃, Ba, Pb, Sr remain retained

No retention of U, Th, Pt, Ir,...

Ra eluted in 3M HNO₃


Further Ba removal via TK102 possible

Tl and Ba eluted in 8M HNO₃




TK102 - Ra purification / recycling

Separation on 1 mL TK102 Resin (100 - 200 μ m) - ~0.5BV/min



Separation on 1 mL SR Resin (100 - 150 μ m) - 0.5BV/min



- SR Resin: high Ba breakthrough starts after 7 – 8 bed volumes
- TK102 Resin: significantly lower Ba breakthrough
- TK102 shows less bleeding than SR Resin




Tc-99m - TK202 Resin

Based on Polyethylene Glycol (PEG)


Less swelling/shrinking than crosslinked PEG

Aqueous Biphasic System (ABS)


Retention of chaotropic anions like TcO_4^- in presence of kosmotrophic anions (SO_4^{2-} , CO_3^{2-} , OH^- , MoO_4^{2-} , ...)



- Separation of Tc-99m from high masses of Mo
- Separation of Re from W (and Ta) possible, too




Dw values for Tc, Re and Mo on TK202 Resin, at varying NaOH concentrations. Tc data taken from Cieszykowska et al.(2).



Dw values for Tc in 5M NaOH using 40 mg TK202 Resin,⁴⁵ increasing amounts of Mo. Data taken from Cieszykowska et



Tc-99m separation from Mo targets – suggested scheme (similar to Zeisler et al.)



- ① Pre-cond. TK202 – 5-7M NaOH → alkaline waste
- ② Load Mo/Tc on TK202 → Mo recovery
- ③ Rinse TK202 – 5-7M NaOH → Mo recovery
- ④ Purge TK202 – Air → Mo recovery
- ⑤ Pre-cond. SCX – HCl then H₂O → Aq. waste
- ⑥ Pre-cond. Al₂O₃ – H₂O → Aq. waste
- ⑦ Elute Tc from TK202 on SCX and load on Al₂O₃ – H₂O
- ⑧ Elute Tc from Al₂O₃ – NaCl 0,9% → Tc recovery

TK202 : 35-75 or 75-150μm
 X : 3-ways valve
 X : 4-ways valve
 SCX : Strong Cation Exchange
 Al₂O₃ : Acidic Alumina

Developed with ReO₄⁻ as TcO₄⁻ surrogate


Re recovered on saline solution from alkaline

Separation with 2g Mo → From 20mL to 2mL

Separation with 200g Mo → From 3L to 20mL



Tc-99m from large Mo targets




Test with 200g Mo

~160 mL TK202 column

Load from 6 - 7M NaOH - elution in water

Pass through C8 cartridge for acidification and Na removal

Final concentration/conversion to 0.9% NaCl or 8 mL AlOxA cartridge





DGA Sheets



TO-DGA (normal DGA) and TEH-DGA (branched DGA) impregnated TLC paper


- Developed at CVUT (Kozempel et al.)

QC of radionuclides and generator eluents


(p.ex. Ra-223, Ac-225/Bi-213, Pb-212, Ge-68/Ga-68 ...)

- TLC scanner or radiometer/LSC or HPGe after cutting


Run under acidic conditions => radionuclidic purity



A scheme of chromatographic separation of mixture of ^{227}Ac and his daughter's nuclides. ^{227}Th remains on start, ^{227}Ac has the retention factor ca 0.2, ^{211}Pb ca 0.7 and ^{223}Ra ca 0.9.



Radiochromatogram measured immediately after separation. Low abundant radiations of ^{227}Ac were not detected.



Radiochromatogram measured one hour after separation. Decay and ingrowth of ^{211}Pb is clearly visible.

More types of sheets under development (selectivities, geometry, support)

- ZR, TK201, ...
- 2D TLC for radionuclide screening ?



CU Sheets

Poster presented at Terachem 2022

(Svedjehed et al.)

QC of Cu radiolabeled peptides (labeled vs free Cu)


- Shown: $[^{61}\text{Cu}]\text{Cu-NOTA-octreotide}$

Spotting/run on three different papers after labeling:

- Whatman and iTLC without modification and
- CU extractant impregnated iTLC paper.

Both iTLC paper (impregnated/non-impregnated) developed in less than 10min, Whatman took 25 – 30 min.

CU extractant impregnated iTLC paper showed superior resolution



- Other systems under development/testing



Some other on-going projects

- Ac rapid purification and 'resalting' via TK222, TK221
- Improvement of radiolysis stability
- Additional 'Sheets'
- Further upscale of radiolanthanide separations
- Scandium separation
 - TK200, TK221, TK222
- Improved Cu separation
- Other radiometals
 - Mn, V, In, Auger (Sb, Pd, Hg, Ag,...)
- At separation
 - TK400, Rn-211/At-211 generator,...
- Decontamination
 - Effluents and reaction wastes
- Fate' of RN in the environment
 - Separation methods
 - Mainly longer lived RN (=> therapy)
 - Ac-225/7, Lu-177(m), radioiodine,...
 - Quantification



Thank you for your attention!



SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Info here

