INHALT

●	Edito	S.1
●	TK222 Resin	S.1-3
●	TK225 Resin	S.4
●	Neuigkeiten	S.5
●	Agenda	S.5
●	Coming soon	S.6

● Edito :

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Zhongtang Wang, Zhaoya Huang, Yun Xie, et al. Method for determination of Pu isotopes in soil and sediment samples by inductively coupled plasma mass spectrometry after simple chemical separation using TK200 resin, *Analytica Chimica Acta*, 1090, 2019, 151-158: www.sciencedirect.com/science/article/pii/S000326701931030X

Zhao Huang, Xiaolin Hou, Xue Zhao, Rapid and Simultaneous Determination of ²³⁸Pu, ²³⁹Pu, ²⁴⁰Pu, and ²⁴¹Pu in Samples with High-Level Uranium Using ICP-MS/MS and Extraction Chromatography, *Anal. Chem.* 2023, 95, 34, 12931–12939: <https://doi.org/10.1021/acs.analchem.3c02526>

Papp, I., Vajda, N. & Happel, S., An improved rapid method for the determination of actinides in water. *J Radioanal Nucl Chem* 331, 3835–3846 (2022): <https://doi.org/10.1007/s10967-022-08389-9>

Ling Zhang, Emilia Vassileva, Determination of ultra-trace level ²⁴¹Am in marine sediment and seawater by combining TK200-TK221 tandem-column extraction chromatography and SF ICP-MS, *Talanta*, 271, 2024, 125724: <https://doi.org/10.1016/j.talanta.2024.125724>

Youyi Ni, Wenting Bu, Ke Xiong, Sheng Hu, Chuting Yang, Liguao Cao. A novel strategy for Pu determination in water samples by automated separation in combination with direct ICP-MS/MS measurement, *Talanta*, 262, 2023, 124710.:

Meet us at our booth F9 at the

Annual Congress of the European Association of Nuclear Medicine (EANM '24)

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- 10th International Conference on Nuclear- and Radiochemistry (NRC10), August 25 – 30, Brighton (GB)
<https://www.rsc.org/events/detail/38385/10th-international-conference-on-nuclear-and-radiochemistry-nrc10>

- 19th International Workshop on Targetry and Target Chemistry (WTTC 19), August 25 – 30, Heidelberg (Deutschland)
https://www.dkfz.de/en/Radiopharmazeutika-Praeklinische-Studien/WTTC19_welcome.html

- Prismap Radiolanthanides Workshop, September 03 – 05, Villigen – PSI (Schweitz) <https://indico.psi.ch/event/15961/page/2872-scope>

- RadWorkshop 2024, September 09 – 13, Riso (Dänemark)
<https://www.conferencemanager.dk/radworkshop2024>

- 3rd Global Meeting COST-NOAR, October 01 – 03, Nantes (Frankreich)
<https://astatine-net.eu/events/cost-third-global-meeting-from-1-to-3-october-2024-in-nantes/>

- “Tumor targeting, Imaging, Radiotherapies” Workshop, October 09 – 12, Erquy (Frankreich) <https://www.cgo-workshop-vecto.fr/>

- 37th Congress of the European Association of Nuclear Medicine, October 19 – 23, Hamburg (Deutschland) <https://eanm24.eanm.org/>

- 67th Annual Radiobioassay and Radiochemical Measurements Conference (RRMC), October 21 – 25, Purdue University, West Lafayette, IN (USA) <https://www.rrmc.co/>

- 6th International Conference on Radioecology & Environmental Radioactivity (ICRER), November 24 – 29, Marseille (Frankreich) <https://www.icrer2024.org/en/>

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


Literatur:

- Baudat, E., Gautier, C., Bagán, H. et al. Optimization of a new radiochemical method based on extraction chromatographic resins and plastic scintillation for measurement of ⁹⁰Sr in nuclear waste. *J Radioanal Nucl Chem.* <https://doi.org/10.1007/s10967-024-09396-8>. 2024
- *I. Giménez, J. Rotger, E. Apellániz, H. Bagán, J. Tent, A. Rigol, A. Tarancón. A new method based on selective fluorescent polymers (PSresin) for the analysis of ⁹⁰Sr in presence of ²¹⁰Pb in environmental samples. *Applied Radiation and Isotopes*, Volume 199, 110879. <https://doi.org/10.1016/j.apradiso.2023.110879>. 2023.
- *Marina Sáez-Muñoz, M.; Bagán, H.; Tarancón, A.; García, J.F.; Ortiz, J.; Carlos, S.; Martorell, S. Rapid methods for radiostrontium determination in aerosol filters and vegetation in emergency situations using PS resin. *Journal of Radioanalytical and Nuclear Chemistry*, 322:1397-1408. <https://doi.org/10.1007/s10967-019-06779-0>. 2019.
- Marina Sáez-Muñoz, M.; Bagán, H.; Tarancón, A.; García, J.F.; Ortiz, J.; Martorell, S. Rapid method for radiostrontium determination in milk in emergency situations using PS resin. *Journal of Radioanalytical and Nuclear Chemistry*. 315, 543–555. 2018
- H. Bagán, A. Tarancón, G. Rauret, J.F. García. Radiostrontium separation and measurement in a single step using plastic scintillators plus selective extractants. Application to aqueous sample analysis. *Analytica Chimica Acta*, 686, 1-2, 50-56. 2011.