



PRODUCT SHEET DISCS

Stainless Steel discs

Main Application:

Sample preparation for alphaspectrometry by electrodeposition

Order N°:

AC-D100-IN25, pack of 100 discs

AC-D250-IN25, pack of 250 discs

AC-D500-IN25, pack of 500 discs

Disc diameter: 25.4 mm

Disc thickness: 0.4 mm

One side of discs protected with vinyl tape

Material: Stainless steel type 304

Composition : Uniloy® 304 type (supplier information) :

Element	Composition (%)
Carbon	0.08 max
Manganese	2.00 max
Phosphor	0.045 max
Sulfur	0.030 max
Silicone	1.00 max
Chrome	18.00-20.00
Nickel	8.00-10.00
Iron	Balance

Silver Disc

Main Application:

Sample preparation for alphaspectrometry by autodeposition

Order N°:

AC-D50-AG25, pack of 25 discs

Disc diameter: 25.4 mm

Disc thickness: 0.4 mm

Both sides of discs protected with vinyl tape

Material: Silver (standard 925)

Composition (supplier information) :

Element	Composition (%)
Silver	92,8
Copper	7,2

Disques Nickel

Main Application:

Sample preparation for alphaspectrometry by autodeposition

Order N°:

AC-D50-IN25, pack of 25 discs

Disc diameter: 25.4 mm

Disc thickness: 0.4 mm

Material: Nickel

Composition (supplier information) :

Element	Composition (%)
Nickel	99,54
Impurities	Balance